

AUSTENAL

AUG 10 2000

4101 W. 51ST STREET
CHICAGO, IL 60632-4287
TEL 773.735.0600
800.621.0381
FAX 773.735.3940

K001799

**DC-TELL for the DCS Precident CAD/CAM System
510(k) Premarket Notification Summary**

The DC-TELL for the DCS Precident CAD/CAM System is a glass fibre reinforced polyamide furnished in a block form that is designed for use as the feed stock for the DCS Precident CAD/CAM System to produce copings and/or substrates for fixed all ceramic dental restorations; i.e. composite bonded to DC-TELL restorations.

DC-TELL has been in use successfully as part of the DCS Precident System in Europe and predicate products exist on the US Market.

The safety and efficacy of DC-TELL has been confirmed by cytotoxicity and allergenic sensitization testing and by experience in with the material in Europe. The suitability of using reinforced polymers for the substructure of dental restorations has also been demonstrated in the dental industry in predicate products.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 10 2000

Mr. Ronald Dudek
Director of Technical Resources
Austenal, Incorporated
4101 West 51st Street
Chicago, Illinois 60632-4287

Re: K001799
Trade Name: Dc-Tell
Regulatory Class: II
Product Code: EBF
Dated: June 7, 2000
Received: June 14, 2000

Dear Mr. Dudek:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Dudek

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Tim S-L/for".

Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): 001799

Device Name: DC-Tell

Indications For Use:

The indications for use for DC-Tell are as substructures for fixed restorations (crowns and bridges) that are veneered with a composite tooth colored material. The substructures are machined using the Precident CAD/CAM System.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)

Gerald W. Smith ^{for MSR}
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K001799