

AUG 9 2000

510(k) Summary
Special 510(k): Device Modification
Stöckert Compact System Version 2.0
(per 21 CFR 807.92)

K002118

1. SPONSOR/APPLICANT

Stöckert Instrumente GmbH
Lilienthalallee 5-7
80939 Munich
Germany

Contact: Helmut Höfl, Director, Quality Assurance and Regulatory Affairs
Telephone: 011 49 89 323 010
Facsimile: 011 49 89 323 4238

Date Prepared: July 12, 2000

2. DEVICE NAME

Proprietary Name: Stöckert Compact System Version 2.0
(SC System V2.0)
Common/Usual Name: Heart-Lung Machine
Classification Names: Multiple (See Table E-1)

Table E-1. SC System V2.0 Classifications

Classification Name	21 CFR	ProCode
Cardiopulmonary bypass heart-lung console	870.4220	74DTQ
Roller type cardiopulmonary bypass blood pump	870.4370	74DWB
Cardiopulmonary bypass pump speed control	870.4380	74DWA
Cardiopulmonary bypass bubble detector	870.4205	74KRL
Cardiopulmonary bypass level sensing monitor/control	870.4340	74DTW
Accessory to the cardiopulmonary bypass console: Temperature Monitor	870.4220	74DTQ
Accessory to the cardiopulmonary bypass console: Timer	870.4220	74DTQ
Cardiopulmonary bypass coronary pressure gauge	870.4310	74DXS

3. PREDICATE DEVICE

Stöckert Compact System, K982014

4. INTENDED USE

The SC System V2.0 is an integrated heart-lung machine consisting of pumps, monitoring, and control elements on a single chassis. It is indicated for speed-controlled pumping through the cardiopulmonary bypass circuit for typical durations of six hours or less, left ventricular venting, cardiotomy suction, administration of cardioplegia solution, when used by a qualified perfusionist who is experienced in the operation of the SC System V2.0.

The SC System has been qualified only for durations appropriate to cardiopulmonary bypass procedures and has not been qualified through in vitro, in vivo, or clinical studies, for long term use as a bridge to transplant, pending recovery of the natural heart, or extracorporeal membrane oxygenation (ECMO) procedures.

5. DEVICE DESCRIPTION

The proposed SC System V2.0 is identical in intended use and fundamental technology to the parent Stöckert Compact System. Modifications are limited to software alterations to the display module of the cardioplegia control unit. The software modifications are as follows:

- Introduction of a control feature to prevent the control module from being switched off, or the flow ratio changed, unless the speed(s) of the connected pump(s) is zero
- Expansion of the selectable flow ratios to include 12/1, 14/1, and 16/1

No other modifications were made to the hardware or software of any Stöckert Compact System component to produce the SC System V2.0.

6. BASIS FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE

The SC System V2.0 is a modification of the Stöckert Compact System and is therefore substantially equivalent to the Stöckert Compact System. This determination is based on equivalence in intended use and technological

characteristics (design and operation). The software modifications described in Section 5 have been validated according to Stöckert Instrumente Design Control procedures, in compliance with the Quality Systems Regulations. Stöckert Instrumente GmbH also believes that any differences between the proposed and parent control modules are minor and raise no new issues of safety or effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

AUG 9 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Medical Device Consultants, Inc.
C/O Ms. Cynthia J. M. Nolte, Ph.D., RAC
Staff Consultant
49 Plain Street
North Attleboro, MA 02760

Re: K002118
Stockert Compact System Version 2.0
Regulatory Class: II (two)
Product Code: DTQ
Dated: July 12, 2000
Received: July 13, 2000

Dear Dr. Nolte:

We have reviewed your Section 510(k) notification of ~~intent to market~~ the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, ~~subject to~~ the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be ~~subject to~~ such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Cynthia J. M. Nolte, Ph.D., RAC

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in dark ink, appearing to read "James E. Dillard III".

James E. Dillard III
Director
Division Of Cardiovascular and
Respiratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K 00 2118

510(k) Number (if known):

Device Name: Stöckert Compact System Version 2.0


Indications For Use:

The Stöckert Compact System Version 2.0 (SC System V2.0) is an integrated heart-lung machine consisting of pumps and monitoring and control elements on a single chassis. It is indicated for speed-controlled pumping through the cardiopulmonary bypass circuit for typical durations of six hours or less, left ventricular venting, cardiomy suction, administration of cardioplegia solution, when used by a qualified perfusionist who is experienced in the operation of the SC System V2.0.

The SC System V2.0 has been qualified only for durations appropriate to cardiopulmonary bypass procedures and has not been qualified through in vitro, in vivo, or clinical studies, for long term use as a bridge to transplant, pending recovery of the natural heart, or extracorporeal membrane oxygenation (ECMO) procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Cardiovascular, Respiratory,
and Neurological Devices

510(k) Number K 00 2118

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐