

*pr/r***5.10 (k) Summary**

August 1, 2000

Contact Person: Bruce Ruefer
 Classification Name: Surgical Mesh
 Common Name: Surgical Mesh
 Trade Name: FluoroTex™ II Surgical Mesh

The FluoroTex II Surgical Mesh is substantially equivalent to Gore-Tex® Soft Tissue Patch. GORE-TEX Soft Tissue Patch and the FluoroTex™ Surgical Mesh consists of a sheet of porous expanded polytetrafluoroethylene (ePTFE); FluoroTex II Surgical Mesh consists of porous expanded polytetrafluoroethylene (ePTFE) reinforced with fluorinated ethylene propylene (FEP). The FluoroTex™ II Surgical Mesh and the predicate devices are intended for the repair of hernias and soft tissue.

Summary of Technological Characteristics


	FluoroTex II Surgical Mesh	GORE-TEX Soft Tissue Patch	FluoroTex Surgical Mesh
Dimensions	1 millimeter and 2 millimeter thickness in a variety of sizes.	1 millimeter and 2 millimeter thickness in a variety of sizes.	1 millimeter and 2 millimeter thickness in a variety of sizes.
Porosity	pore size about 5 to 500 microns	pore size about 10 to 30 microns	pore size about 10 to 30 microns
Material Composition	ePTFE reinforced with FEP	100% ePTFE	ePTFE reinforced with FEP
Material Strength (Kg/cm, 1 millimeter thick material)	16.1 ⁽³⁾	14.8 ⁽¹⁾	13.8 ⁽²⁾
Suture Retention Strength (Kg/pin 1 millimeter thick material)	2.1 ⁽³⁾	1.9 ⁽¹⁾	2.2 ⁽²⁾

Test Notes:

- (1) Reported in Gore literature; n=15; Standard ASTM methods.
 (2) Tested at Bridger Biomed, Inc. labs; n=15; Standard ASTM methods.
 (2) Tested at Bridger Biomed, Inc. labs; n=12; Standard ASTM methods.

5.10 (k) Summary (page 2)**CONCLUSION:**

Mechanical and chemical tests, including material strength and chemical identification of the materials demonstrate that the FluoroTex II Surgical Mesh and the predicate devices, the Gore-Tex Soft Tissue Patch FluoroTex Surgical Mesh, are substantially equivalent.


Bruce G. Ruefer, President

8-1-2000
date

Gore-Tex is a Registered Trademark of W.L. Gore and Associates
 FluoroTex is a Trademark of Bridger Biomed Inc.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

AUG 25 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Bruce G. Ruefer
President
Bridger Biomed, Inc.
2430 N. 7th Street, Suite 4
Bozeman, Montana 59715

Re: K002351
Trade Name: FluoroTex™ II Surgical Mesh
Regulatory Class: II
Product Code: FTL
Dated: August 1, 2000
Received: August 2, 2000

Dear Mr. Ruefer:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Bruce G. Ruefer

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



for Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use:

The FluoroTex™ II Surgical Mesh is intended to be used for the repair of soft tissue and the reconstruction of hernias.

Prescription Use X
(Per 21 CFR 801.109)

Thunell Sayer
(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K002351