

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

NAME OF FIRM: DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
Warsaw, IN 46581-0988

510(K) CONTACT: Arlene C. Saull, RAC
Sr. Regulatory Associate
Phone: 219-371-4904
FAX: 219-371-4940

TRADE NAME: Fenning Hip Stem

COMMON NAME: Femoral Hip Prosthesis

CLASSIFICATION: Class II, per 21 CFR, 888.3358: Hip joint metal/polymer/metal semi-constrained porous-coated uncemented prosthesis.

DEVICE PRODUCT CODE: 87 LPH Prosthesis, Hip, Semi-Constrained, Metal/Polymer, Porous Uncemented

SUBSTANTIALLY EQUIVALENT DEVICES: DePuy AML MMA with Proximal Porocoat
DePuy Titanium Femoral Hip Stem
Stability Femoral Hip Stem

DEVICE DESCRIPTION:

The Fenning Hip Stem is manufactured from titanium alloy. It is a proximally porous-coated, collared, straight-stem design, with distal flutes, a coronal slot, a bullet-shaped distal tip, and is designed with a self-locking taper for use with a DePuy femoral head.

INDICATIONS AND INTENDED USE:

Intended Use:

The Fenning Hip Stem is intended for cementless use with fixation provided by biological tissue ingrowth into the porous coating or for cemented use in which the porous coating serves as a means to augment the fixation of the prosthesis to the bone cement.

Continued on next page.

000004

Indications:

Total hip arthroplasty is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components. Total hip replacement is indicated in the following conditions:

1. A severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia.
2. Avascular necrosis of the femoral head.
3. Acute traumatic fracture of the femoral head or neck.
4. Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemiarthroplasty, surface replacement arthroplasty, or total hip replacement.
5. Certain cases of ankylosis.

BASIS OF SUBSTANTIAL EQUIVALENCE:

DePuy bases its substantial equivalency claim on the design characteristics, which have all been previously cleared by the 510(k) process for one, two, or all three predicate devices. The subject's design includes proximal porous coating, bullet distal tip, distal flutes, coronal slot, modular self-locking taper, straight stem design, modified medial aspect, titanium material, and has a large collar. The subject device and all of the predicates have the same indications, and use the same manufacturing and sterilization methods.

The subject's size range is within the range cleared for the predicate devices.

000005



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

SEP -1 2000

Ms. Arlene C. Saull, RAC
Senior Regulatory Associate
Depuy Orthopaedic, Inc.
700 Orthopaedic Drive
Warsaw, Indiana 46581-0988

Re: K002441
Trade Name: Fenning Hip Stem
Regulatory Class: II
Product Code: LZO and JDI
Dated: August 7, 2000
Received: August 9, 2000

Dear Ms. Saull:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Arlene C. Saull, RAC

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



 Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative and
Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS

510(k) Premarket Notification

510(k) Number (if known) K002441

Device Name Fenning Hip System

Indications for Use:

Total hip arthroplasty is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components. Total hip replacement is indicated in the following conditions:

1. A severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia.
2. Avascular necrosis of the femoral head.
3. Acute traumatic fracture of the femoral head or neck.
4. Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemiarthroplasty, surface replacement arthroplasty, or total hip replacement.
5. Certain cases of ankylosis.

Concurrence of CDRH, Office of Device Evaluation:

Dan R. Lockner
(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K002441

Prescription Use X OR Over-The Counter Use ____ (Per 21 CFR 801.109)

000003