

SUMMARY OF SAFETY AND EFFECTIVENESS

K002678

Sponsor: Biomet, Inc.
56 East Bell Drive
Airport Industrial Park
Warsaw, IN 46582

SEP 19 2000

Contact: Dalene T. Binkley
Phone: (219) 372-1612

Device(s): Ascent™ Knee System:
1. Ascent™ Porous Open Box PS Femoral Component
2. Ascent™ PS Distal Femoral Pegs

Classification: Knee joint patellofemorotibial polymer/metal/polymer semi-constrained cemented prosthesis (CFR 888.3560).

Device Description: There are two components introduced into the Ascent™ Knee System which are as follows: the porous coated open box posterior stabilized (PS) femoral component and the PS distal modular femoral pegs. The indications for use of the Ascent™ Knee System are for the painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved; the correction of varus, valgus, or posttraumatic deformity; and the correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure. The Ascent™ Knee System is for use with bone cement.

Ascent™ Porous Open Box PS Femoral Component

The Ascent™ Porous Open Box PS Femoral Component is the same as the predicate, the Ascent™ Open Box PS Femoral Component, except it has the addition of Ti-6Al-4V plasma spray porous coating to its interior, and it has the addition of a larger size.

The purpose of the porous finish is to enhance cement fixation. The purpose of the additional size is to accommodate larger femurs.

Ascent™ PS Distal Pegs

The Ascent™ PS Distal Femoral Pegs are intended for use with the Ascent™ Open Box PS Femoral Components - porous or Interlok finishes. The modular pegs can fasten into the existing distal augment holes of the components and act as additional alignment devices.

Potential Risks: The potential risks associated with this device are the same as with any joint replacement device. These include, but are not limited to:

| | | |
|-----------------------------|----------------------------|----------------|
| Reaction to bone cement | Blood vessel damage | Bone fracture |
| Deformity of the joint | Soft tissue imbalance | Infection |
| Cardiovascular disease | Delayed wound healing | Hematoma |
| Fracture of the cement | Metal sensitivity | Dislocation |
| Implant loosening/migration | Fracture of the components | Excessive wear |
| Tissue growth failure | Nerve damage | |

000033



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 19 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Dalene T. Binkley
Biomet Orthopedics, Inc
P.O. Box 587
Warsaw, Indiana 46581-0587

Re: K002678
Trade Name: Ascent™ Porous Open Box Posterior Stabilized (PS) Femoral Component
Regulatory Class: II
Product Code: JWH
Dated: August 22, 2000
Received: August 28, 2000

Dear Ms. Binkley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Dalene T. Binkley

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Danne R. Witten

CM Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative and
Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510 (k) NUMBER (IF KNOWN): 4002678

DEVICE NAME: Ascent™ Knee System

INDICATIONS FOR USE:

The indications for use of the Ascent™ Knee are for the painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved; the correction of varus, valgus, or posttraumatic deformity; and the correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure.

This device is for use with bone cement.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED.)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Danner R. Kochner
(Division Sign-Off)

Division of General Restorative Devices

Prescription Use X
(Per 21 CFR 801.109)

OR

510(k) Number 4002678

Over-The-Counter-Use _____
(Optional Format 1-2-96)

000008