

MAY - 1 2001

510(k) Sub. - Ultra Care 2 Natural Rubber Latex Examination Glove  
Made from Allotex™ an enzyme treated natural rubber latex  
with a Protein Content Claim of 200 micrograms or less  
with 15 milligrams or less of Total Particulate per Glove.  
Submission Date: April 2001  
510(k) Number: K002718

[TAB #12]

Attachment #10

Summary of 510(k) Submission

A. INFORMATION

1. SUBMITTER'S  
NAME:

TILLOTSON HEALTHCARE  
CORPORATION

ADDRESS:

360 Route 101  
Bedford, NH 03110 U.S.A.

TELEPHONE NUMBER:

(603) 472-6600

CONTACT PERSON:

F.W. Perrella

DATE SUMMARY PREPARED:

April 2001

2. NAME OF DEVICE  
TRADE OR PROPRIETARY NAME:

Ultra Care 2 Latex Examination Gloves  
made from Allotex™ (enzyme treated)  
Natural Rubber Latex with a Protein  
Content Labeling Claim (200 micrograms or  
less) with 15 milligrams or less of Total  
Particulate per Glove

COMMON OR USUAL NAME:

Examination Glove

CLASSIFICATION  
NAME:

Examination Glove

3. PREDICATE DEVICE IDENTIFICATION  
NAME, NUMBER

1. Ultra Care  
Examination Glove K960247

4. DESCRIPTION OF DEVICE

a. HOW THE DEVICE FUNCTIONS:

Natural Rubber Latex films form a barrier to body fluids and bloodborne pathogens.

b. SCIENTIFIC CONCEPTS THAT FORM THE BASIS FOR THE DEVICE:

The latex rubber is water tight under normal conditions of use. It's tensile properties cause it to conform to the hand, allowing movements necessary for a medical procedure.

c. PHYSICAL AND PERFORMANCE CHARACTERISTICS SUCH AS DESIGN,  
MATERIALS  
AND PHYSICAL PROPERTIES:

Natural Rubber Latex is known to create a barrier to bloodborne pathogens and body fluids. ASTM conforming tensile properties create a glove that is strong

K002718

510(k) Sub. Ultra Care 2 Natural Rubber Latex Examination Glove  
 Made from Allotex™ an enzyme treated natural rubber latex  
 with a Protein Content Claim of 200 micrograms or less  
 with 15 milligrams or less of Total Particulate per Glove.  
 Submission Date: April 2001  
 510(k) Number: K002718

and flexible. The leaching process removes traces of chemical accelerants that may be chemically irritating. The glove is manufactured in accordance with the requirements of ASTM D3578-99 and ASTM D5151-99 requirements.

5. STATEMENT OF INTENDED USE, INCLUDING DESCRIPTION OF THE DISEASES OR CONDITIONS THAT THE DEVICE WILL ADDRESS

This is a disposable device, intended for medical purposes, that is worn on the examiner's hand to prevent contamination between patient and examiner. Examination gloves with protein content labeling are suitable in situations where healthcare worker or patient allergic sensitivity may be a factor.

6. EXPLANATION OF SIMILARITIES OR DIFFERENCES TO PREDICATE DEVICE

- ☐ The modified product has a raw material change whereby the natural rubber latex is treated with proteolytic enzymes to digest natural rubber latex proteins compared to the predicate product.
- ☐ It has no starch donning powder added in the same way as the predicate product with a synthetic inner coating, but with a protein content labeling claim, and Made from Allotex an enzyme treated natural rubber latex claim.

B. IF THE DECISION BASED ON PERFORMANCE DATA

1. DISCUSSION OF NON-CLINICAL TESTS

SPECIFICATION	PROPOSED	PREDICATE
	Synthetic inner coating with no starch donning powder added with Protein Content labeling Claim and Made from Allotex an enzyme treated natural rubber latex claim	Synthetic inner coating with no starch donning powder added
PERFORMANCE STANDARDS	ASTM D3578-99	ASTM D3578-95
WATER TIGHTNESS	ASTM D5151-99	ASTM D5151-92
RESIDUAL PROTEIN	ASTM D5712-99	

2. DISCUSSION OF CLINICAL TESTS

SPECIFICATION	PROPOSED	PREDICATE
<u>SAFETY</u>		
RABBIT IRRITATION	Passes	Passes
GUINEA PIG SENSITIZATION	Passes	Passes

K002718

510(k) Sub. - Ultra Care 2 Natural Rubber Latex Examination Glove  
Made from Allotex™ an enzyme treated natural rubber latex  
with a Protein Content Claim of 200 micrograms or less  
with 15 milligrams or less of Total Particulate per Glove.  
Submission Date: April 2001  
510(k) Number: K002718

3. CONCLUSIONS DRAWN FROM NONCLINICAL AND CLINICAL TESTS THAT  
DEMONSTRATE  
SAFETY EFFECTIVENESS, AND PERFORMANCE => PREDICATE PRODUCT

The Ultra Care 2, Examination Glove has been carefully compared to legally marketed devices in the 510(k). The data summaries indicate that the proposed product meets or exceeds acceptable scores for the predicate product in nonclinical tests, and satisfies the requirements for a safe and effective, no starch donning powder added with 15 milligrams or less of total particulate with protein content labeling claim (200 micrograms or less) per glove and Made from Allotex an enzyme treated natural rubber latex claim medical glove

---

Pursuant to 21 C.F.R. 807.87 (j), I, F.W. Perrella, Ph.D., Vice President R&D certify that to the best of my knowledge and belief and based upon the data and information submitted to me in the course of my responsibilities as the V.P. R&D for TILLOTSON HEALTHCARE CORPORATION, and in reliance thereupon, the data and information submitted in this of the substantial equivalence of this device have been knowingly omitted from this Submission.

---

F.W. Perrella, Ph.D.  
Vice President R&D





DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY - 1 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Frank W. Perrella  
V. P. of Research and Development  
Tillotson Healthcare Corporation  
360 Route 101  
Bedford, New Hampshire 03110-5030

Re: K002718  
Trade/Device Name: Ultra Care Latex Examination Gloves  
made from Allotex (enzyme treated) natural latex with a  
Protein Content Labeling Claim (200 micrograms or less)  
with 15 milligrams or less of Total Particulate per Glove  
Regulation Number: 880.6250  
Regulatory Class: I  
Product Code: Lyy  
Dated: January 31, 2001  
Received: February 1, 2001

Dear Mr. Perrella:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to

comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K002718

510(k) Sub. - Ultra Care 2 Natural Rubber Latex Examination Glove  
Made from Allotex™ an enzyme treated natural rubber latex  
with a Protein Content Claim of 200 micrograms or less  
with 15 milligrams or less of Total Particulate per Glove.  
Submission Date: April 2001  
510(k) Number: K002718

- 3.0 **Indications for Use Statement:** Include the following or equivalent Indications for Use page.  
The information, data and labeling claims in the entire the 510(k) submission must support and agree  
with the Indications for Use statement.

**INDICATIONS FOR USE**

Applicant: Tillotson Healthcare Corporation


5 10(k) Number (if known):\* \_\_\_\_\_

Device Name: Ultra Care 2 Latex Examination Gloves made from Allotex™ (enzyme treated)  
Natural Rubber Latex with a Protein Content Labeling Claim (200 micrograms  
or less) with 15 milligrams or less of Total Particulate per Glove

Indications For Use:

The Ultra Care 2 Examination Glove is "a disposable device intended for medical purposes  
that is worn on the examiner's hand to prevent contamination between patient and examiner."  
(21 CFR 880.6250).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)  
Concurrence of CDRH Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Dental, Infection Control,  
and General Hospital Devices

510(k) Number K002718

Prescription Use \_\_\_\_\_ OR Over-The-Counter \_\_\_\_\_  
Per 21 CFR 801.109  
(Optional Format 1-2-96)

For a new submission, do NOT fill in the 510(k) number blank.