

DEC 15 2000

K002747

510(K) SUMMARY

Reflection® Cross-linked UHMWPE Acetabular Components

In accordance with the Food and Drug Administration Rule to implement provisions of the Safe Medical Devices Act of 1990 and in conformance with 21 CFR 807, this information serves as a Summary of Safety and Effectiveness for the use of the Reflection® Cross-linked UHMWPE Acetabular Components.

Submitter's name:	Smith & Nephew, Inc.
Submitter's address:	1450 Brooks Road Memphis, TN 38116
Submitter's telephone number:	901-399-6487
Contact person:	David Henley
Date summary prepared:	December 8, 2000
Trade or proprietary name:	Reflection® Cross-linked UHMWPE Acetabular Components
Common or usual name:	Polyethylene Acetabular Components
Classification name:	21 CFR 888.3350 Prosthesis, Hip, Semi-Constrained, metal/polymer, Cemented – Class II 21 CFR 888.3358 Prosthesis, Hip, Semi-Constrained, metal/polymer, Uncemented – Class II
Device Product Code and Panel Code:	Orthopedics/87/LPH, JDI

Substantially Equivalent Legally Marketed Devices

- Reflection® Cross-linked UHMWPE Acetabular Components – Smith & Nephew, Inc.
- Inter-Op™ Acetabular System, Durasul™ Acetabular Insert – Sulzer Orthopaedics, Inc.
- Duration® II Acetabular Components-Gas Plasma Sterilization – Howmedica, Inc.
- Osteonics® Crossfire™ Polyethylene Acetabular Components – Osteonics Corp.
- DePuy Duraloc® Acetabular Cup System – DePuy Inc.

Device Description:

The intended use, type of interface, and design features of the Reflection® Cross-linked UHMWPE Acetabular Components are substantially equivalent to the subject identical predicate counterparts fabricated from conventional UHMWPE. Reflection® 10 Mrad cross-linked polyethylene (UHMWPE) acetabular liners are intended to be used only with metal (CoCr) femoral heads.

Device Intended Use:

Total hip components are indicated for individuals undergoing primary and revision surgery where other treatments or devices have failed in rehabilitating hips damaged as a result of trauma, inflammatory joint disease such as rheumatoid arthritis, or noninflammatory degenerative joint disease

(NIDJD) or any of its composite diagnoses such as osteoarthritis; avascular necrosis; traumatic arthritis; slipped capital epiphysis; fused hip; fracture of the pelvis and diastrophic variant; old, remote osteomyelitis with an extended drainage-free period; non-union; femoral neck fracture and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques; femoral osteotomy, or Girdlestone resection; fracture dislocation of the hip; and correction of deformity.

Reflection® Cross-linked UHMWPE Acetabular Components are designed for single use only.

Technological Characteristics:

Reflection® Cross-linked UHMWPE Acetabular Components are substantially equivalent to the subject identical predicate counterparts. The intended use, material, and design features of Reflection® Cross-linked UHMWPE Acetabular Components are substantially equivalent to predicate devices. The modification to the manufacturing process for this polyethylene results in a higher cross-linked material. The safety and effectiveness of this cross-linked polyethylene material in acetabular liner/cup applications, as well as the proposed wear claims, are adequately supported by the substantial equivalence information, materials data, and testing results provided in this premarket notification submission. While Reflection® Cross-linked UHMWPE Acetabular Components are not identical to all of the predicates, any differences that may exist do not significantly affect the safety and effectiveness.

Performance characteristics:

Data indicate that the Reflection® Cross-linked UHMWPE Acetabular Components are substantially equivalent to legally marketed devices.

Wear claims:

The following marketing claims will be made for the Reflection® Cross-linked UHMWPE Acetabular Components:

- No free radicals detectable when using electron spin resonance (ESR) on the final product.
- No increase in oxidation level after accelerated aging compared to conventional, non-irradiated UHMWPE as measured by FTIR.
- No detectable wear as measured gravimetrically in anatomic hip simulator testing. This was compared to the identical predicate counterpart fabricated from conventional, non-irradiated UHMWPE.
- Generation of 72% fewer particles than conventional, non-irradiated UHMWPE in anatomic hip simulator testing.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 15 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. David Henley
Regulatory Affairs Specialist II
Smith & Nephew, Inc.
Orthopaedic Division
1450 Brooks Road
Memphis, Tennessee 38116

Re: K002747
Trade Name: Reflection® Cross-Linked UHMWPE Acetabular Components
Regulatory Class: II
Product Code: LPH and JDI
Dated: December 4, 2000
Received: December 6, 2000

Dear Mr. Henley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 – Mr. David Henley

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a stylized flourish at the end. To the left of the signature is a small, handwritten "for" in cursive.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and
Neurological Devices

Office of Device Evaluation


Center for Devices and
Radiological Health

Enclosure

Indications Statement

Reflection® Cross-linked UHMWPE Acetabular Components

Total hip components are indicated for individuals undergoing primary and revision surgery where other treatments or devices have failed in rehabilitating hips damaged as a result of trauma, inflammatory joint disease such as rheumatoid arthritis, or noninflammatory degenerative joint disease (NIDJD) or any of its composite diagnoses such as osteoarthritis; avascular necrosis; traumatic arthritis; slipped capital epiphysis; fused hip; fracture of the pelvis and diastrophic variant; old, remote osteomyelitis with an extended drainage-free period; non-union; femoral neck fracture and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques; femoral osteotomy, or Girdlestone resection; fracture dislocation of the hip; and correction of deformity.


(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K002747