

K003181

NOV - 7 2000

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

**The Implex Continuum Hip System – Acetabular Components**

**Submitter Name:** Implex Corp.  
**Submitter Address:** 80 Commerce Drive  
Allendale, New Jersey 07401-1600  
**Contact Person:** John Schalago  
**Phone Number:** (201) 818-1800  
**Fax Number:** (201) 995-9722  
**Date Prepared:** November 6, 2000  
**Device Trade Name:** The Implex Continuum Hip System  
**Device Common Name:** Acetabular Components  
**Classification Number and Name:** 21 CFR § 888 - various

---

**Substantial  
Equivalence:**

The term "substantial equivalence" as used in this 510(k) notification is limited to the definition of substantial equivalence found in the Federal Food, Drug and Cosmetic Act, as amended and as applied under 21 CFR 807, Subpart E under which a device can be marketed without premarket approval or reclassification. A determination of substantial equivalency under this notification is not intended to have any bearing whatsoever on the resolution of patent infringement suits or any other patent matters. No statements related to, or in support of substantial equivalence herein shall be construed as an admission against interest under the US Patent Laws or their application by the courts.

**Device Description:**

Implex Continuum Hip System – Acetabular Components are offered in primary and revision options. The system includes monoblock, modular, and non-modular acetabular components. These components are manufactured from Hydrocel Porous Tantalum, CP titanium, Titanium alloy, and UHMWPE. The acetabular components are offered in OD sizes 40mm to 70mm and four ID sizes 22, 26, 28 and 32mm. Implex Acetabular components with cross-linked polyethylene are cleared for use with Implex and Zimmer metallic femoral heads only.

The Continuum Hip System Acetabular Components are intended for use with Implex and Zimmer femoral heads and bone screws.

### **510(k) Summary (Continued)**

<b>Indications for Use:</b>	The Implex Continuum Hip System is intended for use where severe degeneration, trauma, or other pathology of the hip joint indicates cemented, cementless or hybrid total hip arthroplasty.
<b>Device Technological Characteristics and Comparison to Predicate Device:</b>	The device technological characteristics are not affected by interfacing of Implex Continuum Hip System - Acetabular components with Zimmer femoral components or the labeling change described herein. The subject device and the predicate devices are manufactured from same materials, have the same indications for use, and utilize same surgical techniques and instrumentation. Implex acetabular components with cross-linked polyethylene are cleared for use with Zimmer/Implex metallic femoral heads only.
<b>Performance Data:</b>	Evaluation of tolerances, tolerance stack-up, and range of motion supports the interface of Implex/Zimmer femoral head or bone screws. The performance data provided is not intended to support the interfacing of Implex acetabular components with cross-linked polyethylene and Zimmer/Implex ceramic bearing heads.
<b>Conclusion:</b>	The Implex Continuum Hip System is substantially equivalent to the identified predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

NOV - 7 2000

Mr. John A. Schalago, RAC  
Manager, Regulatory Affairs  
Implex Corporation  
80 Commerce Drive  
Allendale, New Jersey 07401-1600

Re: K003181

Trade Name: Implex Continuum Hip System - Acetabular Components  
Regulatory Class: II  
Product Code: LPH  
Dated: October 10, 2000  
Received: October 11, 2000

Dear Mr. Schalago, RAC:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Mr. John A. Schalago, RAC

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Celia M. Witten". The signature is written in a cursive, flowing style.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if  
known):K003181

Device Name:

The Implex Continuum Hip System-Acetabular  
Components

Indications For Use:

The Implex Continuum Hip System-Acetabular Components are intended for use where severe degeneration, trauma, or other pathology of the hip joint indicates cemented, cementless or hybrid total hip arthroplasty.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Concurrence of CDRH; Office of Device Evaluation (ODE)

Prescription  
in Use  
(Per 21 CFR 801.109)

OR...

Over-The-  
Counter Use

(Optional Format 1-2-96)

*for Mark N. Mellgren*  
(Division Sign-Off)  
Division of General Restorative Devices  
510(k) Number K003181