

DEC 28 2000

K003212

510(K) Summary

Disc-O-Tech Medical Technologies, Ltd.
Fixion Intramedullary Nailing System – 6.7/10.0 Nail

Company Name

Disc-O-Tech Medical Technologies, Ltd.
3 Hasadnaot St. Herzelia
Israel, 46728

Submitter's Name and Contact Person

Jonathan S. Kahan, Esq.
555 Thirteenth Street, NW
Washington, DC 20004
Tel: 202-637-5794
Fax: 202-637-5910

Date Prepared

August, 2000

Trade/Proprietary Name

FixionTM Intramedullary Nailing System (Fixion IM Nail)

Classification Name

Intramedullary Fixation Rod
21 CFR § 888.3020
Class II

Predicate Devices

1. Fixion Intramedullary Nailing System (K990717) by Disc-O-Tech.
2. True/Flex Upper Extremity IM Nail (K902264) by Encore (Applied Osteo Systems, Inc.)
3. Seidel Humeral Locking Nail (K883882, K924004, K925544, K931256) by Howmedica
4. AO Unreamed Humeral (UHN) Nail (K933518) by Synthes.

Performance Standards

The following standards were used:

1. The Fixion IM Nail is manufactured from 316L Stainless Steel, which meets the requirements of ASTM F138 - Standard Specification for Stainless steel Bar and Wire for Surgical Implants.
2. The 4 point bending mechanical testing was performed according to ASTM F1264 - Standard for Mechanical Performance Considerations for Intramedullary Fixation Devices.
3. The torsional mechanical testing was performed according to ASTM F383 - Standard Practice for Static Bend and Torsion Testing of Intramedullary Rods.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 28 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Johnathan S. Kahan, Esq.
Disc-O-Tech Medical Technologies, LTD.
c/o Hogan & Hartson L.L.P.
Columbia Square
555 Thirteenth Street, N.W.
Washington, DC 20004

Re: K003212

Trade Name: Fixion Intramedullary Nailing System ("Fixion IM Nail")
Regulatory Class: II
Product Code: JDS
Dated: October 13, 2000
Received: October 13, 2000

Dear Mr. Kahan:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Mr. Johnathan S. Kahan, Esq.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for Mark N. Milken

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and
Neurological Devices

Office of Device Evaluation

Center for Devices and
Radiological Health

Enclosure

Indication for Use

510(K) Number (if known): _____

Device Name: Fixion™ Intramedullary Nailing System (Fixion™ IM Nail) – Nail size
6.7/10.0

Indication for Use: Nail 6.7/10.0 of the Fixion™ Intramedullary Nailing System family of nails is intended for use in the fixation of humeral fractures, including diaphyseal fractures. It is indicated for use in shaft fractures 5cm below the surgical neck to 5cm proximal to the distal end of the medullary canal.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(per 21 CFR 801.109)

OR

Over the Counter Use _____

for Mark H. Mulhens
(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K003212