

NOV - 6 2000

K003296

## SUMMARY OF SAFETY AND EFFECTIVENESS

**Sponsor:** Biomet Inc.  
Airport Industrial Park  
P.O. Box 587  
Warsaw, IN 46581-0587

**Contact Person:** Tracy J. Bickel  
(219) 372-1761

**Device(s):** Biomet® SCK Knee System

**Classification Name:**

Knee joint patellofemorotibial polymer/metal/polymer semi-constrained cemented prosthesis (CFR 888.3560)

**Device Classification:** Class II

**Device Product Code:** 87 JWH 888.3560

**Intended Use:**

- Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, and/or traumatic arthritis where one or more compartments are involved.
- Correction of varus, valgus, or posttraumatic deformity.
- Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure

This device is for use with bone cement

**Device Description:**

The Biomet® SCK Knee System has an open box femoral component with locked tibial bearing, that will utilize the Maximum Congruent Knee (K915132) tibial trays. All SCK implants are manufactured from the same materials as the Maximum Congruent Knee (K915132) implants

The tibial bearing has a high, wide post with a lip at the top. The sagittal tibial profile has a 'dished' surface that will interact with the distal radius of the femoral component. The tibial bearings are available in sizes 59, 63/67, 71/75, 79/83, and 87/91, with thickness ranging from 10mm thick to 24mm in 2mm increments. The SCK tibial bearing is composed of ArCOM® UHMWPE (ASTM F648).

The SCK femoral component is made from cast Cobalt Chrome alloy (ASTM F75). The femoral component has an open box structure with an outer width of inches. For stem attachment there is a boss angle of 6° located superiorly to the box. The SCK femoral components are available in five anatomical sizes, 60-80mm in five millimeter increments.

00095

**Potential Risks:** The potential risks associated with this device are the same as with any joint replacement device. These include, but are not limited to:

|                             |                            |                |
|-----------------------------|----------------------------|----------------|
| Reaction to bone cement     | Blood vessel damage        | Bone fracture  |
| Deformity of the joint      | Soft tissue imbalance      | Infection      |
| Cardiovascular disease      | Delayed wound healing      | Hematoma       |
| Fracture of the cement      | Metal sensitivity          | Dislocation    |
| Implant loosening/migration | Fracture of the components | Excessive wear |
| Tissue growth failure       | Nerve damage               |                |

**Predicate Device(s):**

Maximum Congruent Knee, K915132- Manufactured by Biomet, Inc. Warsaw, IN.

00096



NOV - 6 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Tracy J. Bickel  
Regulatory Specialist  
Biomet, Inc.  
P.O. Box 587  
Warsaw, Indiana 46581

Re: K003296  
Trade Name: Biomet SCK Knee System  
Regulatory Class: II  
Product Code: JWH  
Dated: October 18, 2000  
Received: October 20, 2000

Dear Ms. Bickel:

We have reviewed your Section 510(k) notification of intent to market the ~~device referenced~~ above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 – Ms. Tracy J. Bickel

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a stylized flourish at the end.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and  
Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K00-3296  
Device Name: **Biomet® SCK Knee System**  
Indications for Use:

The Biomet® SCK Knee System is indicated for:

- 1) Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, and/or traumatic arthritis where one or more compartments are involved.
- 2) Correction of varus, valgus, or posttraumatic deformity.
- 3) Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure.

This device is for use with bone cement

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use \_\_\_\_\_  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)

*f* Mark N. Miller  
(Division Sign-Off)

Division of General Restorative Devices

510(k) Number K 00 3296

00007