



DEC 21 2000

**Premarket Notification [510(k)] Summary
for Hem-O-Lok® Ligation Clips**

1. Submitter Name, Address, and Date of Submission

Brian Young
Sr. Regulatory Affairs Manager
Weck Closure Systems
One Weck Drive
Research Triangle Park, NC 27709

Telephone: (919) 361-4041

Facsimile: (919) 361-3914

Submitted: October 20, 2000

2. Name of the Device, Common, Proprietary (if known), and Classification

Classification Name: Implantable clip

Common Name: Ligating clip

Proprietary Name: Hem-O-Lok® Ligating clip

3. Identification of the legally marketed device to which the submitter claims equivalence

The Weck Closure Systems Hem-O-Lok® clip is substantially equivalent to Hem-O-Lok® clips cleared under Weck's previous 510(k) filing number K993157, and to preamendment Weck Hemoclip® surgical occluding clips.

4. Description of the Device

Weck Closure System's Hem-O-Lok® ligation clip is a manually applied hemostatic clip intended to connect internal tissues to aid healing. Hem-O-Lok® causes hemostasis through vessel ligation. The clip is nonabsorbable and is manufactured from polyacetal.

The clips are housed in a cartridge and packaged in a rigid plastic blister with Tyvek coated lidding which is sold sterile. The method of sterilization is EtO with a SAL of 10^{-6} . The blister packs are fitted into an overpack carton which serves as the sales unit.

5. Intended Use of the Device

Indications for use: "Hem-O-Lok® ligating clips are intended for use in procedures involving ligation of vessels or tissue structures. Surgeons should apply the appropriate size clip for the size of the vessel or tissue structure to be ligated such that the clip completely encompasses the vessel or tissue structure."

6. Summary of Technological Characteristics

The technological characteristics are the same as or equivalent to the predicate devices. The polyacetal material used in the clips is shown to be biocompatible according to ISO 10993-1 requirements.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 21 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Brian J. Young
Senior Regulatory Affairs Manager
Weck Closure Systems
One Weck Drive
P.O. Box 12600
Research Triangle Park, North Carolina 27709

Re: K003337
Trade Name: Hem-O-Lok® Ligating Clips
Regulatory Class: II
Product Code: FZP
Dated: October 24, 2000
Received: October 25, 2000

Dear Mr. Young:

We have reviewed your Section 510(k) notification of intent to market ~~the device referenced~~ above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Brian J. Young

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,


Celia M. Witten, Ph.D., M.D. *for*
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K003337

2 Statement of indications for use

510(k) Number (if assigned):

Device Name:

Hem-O-Lok® ligating clips

INDICATIONS FOR USE

Hem-O-Lok® ligating clips are intended for use in procedures involving ligation of vessels or tissue structures. Surgeons should apply the appropriate size clip for the size of the vessel or tissue structure to be ligated such that the clip completely encompasses the vessel or tissue structure.

Concurrence of CDRH, Office of Device Evaluation (ODE)

NPD for cmw
(Division Sign-Off)

Division of General Restorative Devices

510(k) Number

K003337

Prescription Use ✓
(Per 21 CFR 801.109)