



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR 15 2001

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Dr. H. Anderson  
Regulatory Affairs Officer  
Randox Laboratories Ltd.,  
Biochemical Manufacturers  
Ardmore, Diamond Road  
Crumlin, Co.  
Antrim,  
United Kingdom  
BT29 4QY

Re: K003346  
Trade Name: Glucose (GOD – PAP)  
Regulatory Class: II  
Product Code: CGA  
Dated: February 14, 2001  
Received: February 16, 2001

Dear Dr. Anderson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical Laboratory Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

510(k) Number (if known) NOT KNOWN K 003346Device Name GLUCOSE (GOD - PAP)LIQUID REAGENT**Indications For Use :**

The Randox Laboratories Ltd Glucose (GOD - PAP) test kit is an *in vitro* diagnostic reagent for the quantitative determination of glucose in serum. The method is based on the enzymatic oxidation of glucose by glucose oxidase. The hydrogen peroxide produced reacts, catalysed by peroxidase, with phenol and 4-aminophenazone to form a red-violet quinoneimine dye. The intensity of the final colour is directly proportional to the glucose concentration and is measured at 505nm. The measurement of glucose in serum is important in the diagnosis and treatment of carbohydrate metabolism disorders such as diabetes mellitus, neonatal hypoglycaemia, and idiopathic hypoglycemia, and of pancreatic islet cell carcinoma.

This Application Sheet has been developed for the Hitachi 717 analyser and must only be used by suitably qualified personnel under appropriate laboratory conditions.

[Signature]  
 (Division Sign-Off)  
 Division of Clinical Laboratory Sciences  
 510(k) Number K 003346

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
 (Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
 (Optional format 1-2-96)