

NOV 30 2000

K003379

Special 510(k) Summary - Device Modification
Summary of Safety and Effectiveness for the
V-40™/C-Taper Adapter Sleeve

Submission Information

Name and Address of the Sponsor
Of the 510(k) Submission

Howmedica Osteonics Corp
59 Route 17
Allendale, NJ 07401-1677

Contact Person:

Karen Ariemma
Regulatory Affairs Specialist

Date of Summary Preparation:

October 19, 2000

Device Identification:

Proprietary Name:

V-40™/C-Taper Adapter Sleeve

Common Name:

Adapter Sleeve for Femoral Head

Classification Name and Reference

Hip Joint, Metal/Ceramic/Polymer, Semi-
Constrained, Cemented or Nonporous
Uncemented Prosthesis, 21 CFR §888.3353

Predicate Device Identification

The V-40™/C-Taper Adapter Sleeve is substantially equivalent to the Osteonics® Titanium Adapter Sleeve which was determined substantially equivalent via 510(k) 885102.

Device Description

The Osteonics® Titanium Adapter Sleeve is a tapered sleeve component with a female Morse taper to provide locking with a Howmedica Osteonics' femoral stem with a Morse taper. In addition, the sleeve has a tapered male exterior surface that provides locking with an Osteonics® C-Taper Alumina Ceramic Head. The V40™/C-Taper Adapter Sleeve is a modification of the predicate Osteonics® Titanium Adapter Sleeve. The modification involves changing the inner taper diameter from a Morse taper to a V40™ taper and removing the cap of the sleeve. This modification is designed to allow both Osteonics® Alumina C-Taper Heads and Zirconia C-Taper Ceramic Heads to mate with Howmedica Osteonics' femoral stems with a V40™ taper. The modified sleeve is identical to the unmodified sleeve in every aspect with the exception of the female taper angle. The subject and predicate devices are fabricated from the same material.

Intended Use

The intended use of the V40™/C-Taper Adapter Sleeve is identical to that of the predicate Osteonics® Titanium Adapter Sleeve. The V40™/C-Taper Adapter Sleeve is intended to allow either an Osteonics® C-Taper Alumina Head or an Osteonics® C-Taper Zirconia Ceramic Head to mate with any Howmedica Osteonics' femoral stem with a V40™ taper. The V40™/C-Taper Adapter Sleeve is a single-use device.

Indications for Use**For Use as a Total Hip Replacement:**

- Painful disabling joint disease of the hip resulting from: degenerative arthritis, rheumatoid arthritis, post-traumatic arthritis or late stage vascular necrosis.
- Revision of previous cup arthroplasty or other procedures
- Clinical management problems where arthrodesis or alternative reconstructive techniques are less likely to achieve satisfactory results.

For Use as a Bipolar Hip Replacement

- Femoral head/neck fractures or non-unions.
- Aseptic necrosis of the femoral head.
- Osteo-, rheumatoid, and post traumatic arthritis of the hip with minimal acetabular involvement or distortion.

Other Considerations:

- Pathological considerations or age considerations which indicate a more conservative acetabular procedure and an avoidance of the use of bone cement in the acetabulum.
- Salvage of failed total hip arthroplasty

Statement of Technological Comparison

The subject V40™/C-Taper Adapter Sleeve share the same material, intended use, and basic design concepts as the predicate Osteonics® Titanium Adapter Sleeve.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 30 2000

Ms. Elizabeth A. Staub
Vice President, Quality Assurance/ Regulatory Compliance/ Clinical Research
Stryker Howmedica Osteonics Corporation
59 Route 17
Allendale, New Jersey 07401-1677

Re: K003379
Trade Name: V40™/C-Taper Adapter Sleeve
Regulatory Class: II
Product Code: LZO
Dated: October 30, 2000
Received: October 31, 2000

Dear Ms. Staub:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Elizabeth A. Staub

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for Mark N. Milken

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and
Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K003379

Device Name: V40™/C-Taper Adapter Sleeve

Indications For Use:

The indications for use of the subject device, in keeping with those of other legally marketed Howmedica Osteonics' adapter sleeves is as follows:

For Use as a Total Hip Replacement:

- Painful disabling joint disease of the hip resulting from: degenerative arthritis, rheumatoid arthritis, post-traumatic arthritis or late stage vascular necrosis.
- Revision of previous cup arthroplasty or other procedures
- Clinical management problems where arthrodesis or alternative reconstructive techniques are less likely to achieve satisfactory results.

For Use as a Bipolar Hip Replacement


- Femoral head/neck fractures or non-unions.
- Aseptic necrosis of the femoral head.
- Osteo-, rheumatoid, and post traumatic arthritis of the hip with minimal acetabular involvement or distortion.

Other Considerations:

- Pathological considerations or age considerations which indicate a more conservative acetabular procedure and an avoidance of the use of bone cement in the acetabulum.
- Salvage of failed total hip arthroplasty

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

for 
(Division Sign-Off)

Division of General Restorative Devices

510(k) Number K003379

Prescription Use X

OR

Over-The-Counter Use

(Per 21 CFR 801.109)

(Optional Format 1-2-96)