

JAN 29 2001

K 003397
10F2

510(k) Summary

Name of Sponsor:

DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
Warsaw, Indiana 46581-0988
Est. Reg. No. 1818910

510(k) Contact:

Marcia J. Arentz
Senior Regulatory Associate
Phone: (219) 371-4944
FAX: (219) 371-4940

Trade Name:

Colles C Series Frame Sterile Pack

Common Name:

External Fixation device

Classification:

Class II Device per 21 CFR 888.3030:
Single/Multiple component metallic bone
fixation appliances and accessories.

Device Product Code:

Code: **87KTT** Appliance, Fixation,
Nail/Blade/Plate combination, multiple
component.
No performance standards have been established
under Section 514 of the Federal Food, Drug,
and Cosmetic Act for femoral hip stems.

Substantially Equivalent Device:

| | |
|--|----------------|
| Threaded Pins Self-Drilling | K781251 |
| Colles Fixator | K842768 |
| Mini-Fixator for Colles Fractures | K982982 |

Device Descriptions:

The Colles C Series Frame Sterile Pack is an external fixation device used in the treatment of fractures of the wrist. The system is comprised of half pins that are implanted through the skin attached to adjustable connecting rods that can apply traction. The kit also includes instrumentation required for the surgery.

000005

510(k) Summary (continued)**Intended use:**

For use in the fixation of distal fractures of the upper extremity. External fixation systems allow control of bone segments including angulation, rotation and displacement.

Indications for use:

Commuted, intra-articular distal radius fractures (Frykman-Classification III-VIII), Bilateral Colles' fracture, failed closed reduction with casting.

Substantial equivalence:

The Colles Fixator is the same device cleared in K842768. The pins are the same design and intended use as those cleared in K781251. The fixator and the pins are packaged together with instruments in a sterile package and are substantially equivalent to the Biomet Mini-Fixator for Colles Fractures cleared in K982982



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JAN 29 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Marcia J. Arentz
Senior Regulatory Associate
DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
P.O. Box 988
Warsaw, Indiana 46581-0988

Re: K003397
Trade Name: Colles C Series Frame Sterile Pack
Regulatory Class: II
Product Code: KTT
Dated: October 30, 2000
Received: November 1, 2000

Dear Ms. Arentz:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for Celia M. Witten, Ph.D., M.D.

Director
Division of General, Restorative and
Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K003397

Device Name: Colles Fixator

Indications for Use:

For use in the fixation of distal fractures of the upper extremity. External fixation systems allow control of bone segments including angulation, rotation and displacement.

Concurrence of CDRH, Office of Device Evaluation

Prescription Use X
(Per 21 CFR 801.109)

OR Over-The-Counter Use

for Mark A. Miller
(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K003397

000003