

FEB 14 2001

K003591  
PAGE 10P

SECTION 5  
SUMMARY OF SAFETY AND EFFECTIVENESS

**Date Summary Prepared:**

19<sup>th</sup> November 2000

**Name of Device:**

Proprietary name:	TensCare XL-Y2
Common name:	TENS device
Classification name:	Stimulator, Nerve, Transcutaneous, for Pain Relief 84GZJ; 21 CFR 882.5890.
Device Classification:	Class II
Predicate Device:	FUJI TENS 804SIII (K893874 B)
Device Description:	A portable TENS device for pain relief.
Intended Purpose/Use:	TENS is used for the relief and management of symptomatic intractable pain and/or as an adjunctive treatment in the management of post-surgical and post traumatic acute pain.
Technological Comparison:	The TENS CARE has basic technological characteristics that are substantially equivalent to the predicate device. The differences in technological characteristics are the use of a Microprocessor for the control of all functions, the use of pre-set output energy levels selectable by depression of a Button (as opposed to rotational control knobs on the predicate device) and the use of 'shrouded patient cable connectors' to comply with FDA's Final Rule "Medical Devices: Establishment of Performance Standards for Electrode Lead Wires and Patient Cables".

K003591  
PAGE 2 OF 2

## SECTION 5 SUMMARY OF SAFETY AND EFFECTIVENESS

Labelling Comparison: The Labelling is substantially equivalent to that of the predicate device.

Non-Clinical Testing: The results of Bench Testing demonstrate that the output characteristics of the TENS CARE XL-Y2 are substantially equivalent to those of the predicate device.

Clinical Testing: Clinical Testing was not necessary as no new or innovative aspects have been introduced.

Further safety information: The "TensCare XL-Y2" device has been on the European Market for the past two years. During this time a review of Customer Complaints, Returned Product and the results of Post Market Feedback, has demonstrated that the product has performed as Intended, to it's Specified Requirements. The data analysed is summarised in this submission and the full data is available upon request. The Certificate of authority to CE Mark the "TensCare XL-Y2" in accordance with the Medical Device Directive 93/42/EEC is included in Section 12 of this submission.

**Conclusions:** The TensCare XL-Y2 is substantially equivalent to the predicate device and any differences between the devices do not pose any new questions of safety and effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB 14 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

TensCare Ltd.  
C/O Bernard Tremaine  
Medical Device & QA Consultancy  
76, Stockport Road  
Timperley  
Cheshire.  
WA15 7SN. United Kingdom

Re: K003591  
Trade Name: TensCare, Model XL-Y2  
Regulatory Class: II  
Product Code: GZJ  
Dated: November 19, 2000  
Received: November 21, 2000

Dear Mr. Tremaine:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 – Mr. Bernard Tremaine

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Mark N. Melkerson", is written over the typed name of Celia M. Witten.

Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): \_\_\_\_\_

Device Name: TENS XL-Y2

Indications For Use:

"TENS XL-Y2" TENS UNIT

"For the symptomatic relief of chronic intractable pain"

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*for Mark N. Miller*  
\_\_\_\_\_  
(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

Prescription Use \_\_\_\_\_  
(Per 21 CFR 801.109)

OR

510(k) Number K003591  
Over-The-Counter Use \_\_\_\_\_