

DEC 20 2000

510(k) Summary of Safety and Effectiveness

ArthroCare Corporation
Visage® Cosmetic Surgery System

K003624

General Information

Manufacturer:

ArthroCare, Corporation
595 North Pastoria Avenue
Sunnyvale, CA 94085-2936

Establishment Registration Number:

2951580

Contact Person:

Betty M. Johnson
Manager, Regulatory Affairs

Date Prepared:

November 16, 2000

Device Description

Classification Name:

Electrosurgical Cutting and Coagulation
Device and Accessories (21 CFR 878.4400)

Trade Name:

Visage® Cosmetic Surgery System

Generic/Common Name:

Electrosurgical Device and Accessories

Predicate Devices

Visage Cosmetic Surgery System

K992180

Intended Uses

The Visage Cosmetic Surgery System is a bipolar electrosurgical device intended for general dermatologic surgery that may include skin resurfacing for the treatment of wrinkles, rhytids, and furrows, as well as soft tissue resection/removal and hemostasis/coagulation. It is intended to be used in procedures using conductive solutions such as normal saline

Product Description

The Visage Cosmetic Surgery System is a bipolar, high frequency electrosurgical system consisting of three components: an electrosurgical generator called the Controller, a disposable Stylet, a reusable Handpiece and Cable Unit, and a reusable Cable Adapter.

Substantial Equivalence

This Special 510(k) proposes a modification in the performance specifications, dimensional specifications, and labeling for the Controller and Handpiece and Cable Unit components of the Visage Cosmetic Surgery System, which was previously cleared in K992180 on March 20, 2000. The technology, principle of operation, and the intended use of the Visage System remains the same as in the predicate cleared 510(k).

Summary of Safety and Effectiveness

The proposed modifications to the Controller and Handpiece and Cable components of the Visage Cosmetic Surgery System, as described in this submission, are substantially equivalent to the predicate device. The proposed modification in performance specifications, labeling, dimensional, and sterilization specifications are not substantial changes or modifications, and do not significantly affect the safety or efficacy of the devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

DEC 20 2000

Ms. Betty Johnson
Manager, Regulatory Affairs
ArthroCare Corporation
595 North Pastoria Avenue
Sunnyvale, California 94085

Re: K003624
Trade Name: Visage® Cosmetic Surgery System
Regulatory Class: II
Product Code: GEI
Dated: November 21, 2000
Received: November 24, 2000

Dear Ms. Johnson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (~~for the indications for~~ use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, ~~or to devices that~~ have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Betty Johnson

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications Statement

Device Name: Visage® Cosmetic Surgery System

510(k) Number: K 003624

Indications for use:

The Visage Cosmetic Surgery System is a bipolar electrosurgical device intended for general dermatologic surgery that may include skin resurfacing for the treatment of wrinkles, rhytids, and furrows, as well as soft tissue resection/removal and hemostasis/coagulation. It is intended to be used in procedures using conductive solutions such as normal saline.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

MM for YMM

(Division Sign-Off)

Division of General Restorative Devices

510(k) Number 003624

Prescription Use

X

OR

Over-the-Counter
Use

(Per 21 CFR
801.109)