

FEB 20 2001

K003659

**510(k) Summary**  
**Bionx Implants Inc.**  
**PLGA Pin**

**Submitter's Name, Address, Telephone Number, and Contact Person**

Bionx Implants, Inc.  
1777 Sentry Parkway West  
Gwynedd Hall, Suite 400  
Bluebell, PA 19422

Contacts: Gerard S. Carlozzi  
President and Chief Executive Officer  
Phone: (215) 643-5000  
Facsimile: (215) 653-0984

Bionx Implants Ltd.  
Tuija Annala  
Quality Manager  
P.O.Box 3  
FIN-33721 Tampere  
Finland  
Phone: 358-3-316 5679  
Facsimile: 358-3-316 5629

**Date prepared:** November 7<sup>th</sup>, 2000

**Name of the device:**

A. Trade or Proprietary Name: PLGA Pin

B	Common Name:	Bioabsorbable bone fixation pin
C.	Classification Name:	Bone Fixation Pin, class II, HTY

### **Predicate Devices:**

The predicate devices are Bionx Implants Inc. Biofix SR-PGA Pin (K890902), SmartPin™ (K925098) and SmartNail™ (K993074). Also our competitors have comparable predicate devices in the market, like Johnson & Johnson Orthopedics, Inc. Orthosorb (K864912, K882979 and K901256) and Synthes (U.S.A) Polypin (K961608).

### **Intended Use:**

PLGA Pin is indicated for fixation of fragments of fractured non-load bearing bones, osteotomies and arthrodeses, for example in the fixation of apical fragments, osteochondral fragments and cancellous/non-load bearing fragments in the presence of appropriate immobilization.

### **Device Description:**

The device description of PLGA Pin is as follows.

- Composed of poly-L/G-lactide copolymer
- Lengths 20, 40 and 60mm
- Diameters 1.1, 1.5, 2.0 and 3.2mm

The dimensions and shape are completely identical with the Biofix SR-PGA Pin (K890902) and SmartPin™ (K925098).

**Substantial Equivalence:**

PLGA Pin has the following similarities to the cleared Bionx Implants Inc. Biofix SR-PGA Pin (K890902), SmartNail (K993074) and SmartPin™ (K925098):

- has the same or similar indicated use
- uses the same operating principle
- incorporates the same basic design
- is manufactured by same machinery
- is packaged and sterilized using the same materials and processes
- has same 2 years shelf life than Biofix SR-PGA Pin (K890902)

PLGA Pins have the following similarities to the cleared Johnson & Johnson Orthopedics, Inc. Orthosorb (K864912, K882979 and K901256) and Synthes (U.S.A) Polypin (K961608):

- has the same or similar indicated use
- use the same operating principle
- incorporate the same basic design

In summary, PLGA Pin described is substantially equivalent to the predicate device. This change of raw material does not raise any problems concerning safety or efficacy of the device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB 20 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Tuija Annala  
Quality Manager  
Bionx Implants, Inc.  
c/o Biounx Implants, LTD  
P.O. Box 3, Fin-33721  
Hermiankatu 6-8 L  
Tampere,  
Finland

Re: K003659  
Trade Name: PLGA Pin  
Regulatory Class: II  
Product Code: HTY  
Dated: November 7, 2000  
Received: November 28, 2000

Dear Ms. Annala:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Tuija Annala

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for Mark A. Melber

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## INDICATIONS FOR USE

510(K) Number (if known): K003659

Device Name: PLGA Pin

### Indications for Use:

PLGA Pin is intended for use in the fixation of fragments of fractured non-load bearing bones, osteotomies and arthrodeses, for example in the fixation of apical fragments, osteochondral fragments and cancellous/non-load bearing fragments.

PLGA Pin is not intended for use in and is contraindicated for: 1) Fractures and osteotomies of cortical bone (except cortical bones of the foot and the hand), 2) Fractures and osteotomies in weight bearing cancellous bone, 3) Situations where internal fixation is otherwise contraindicated, e.g., active or potential infection and where patient co-operation cannot be guaranteed (e.g. alcoholism).

(Please do not write below this line – continue on another page is needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use yes  
(Per 21 CFR 801.109)

OR Over-The-Counter Use No

for Mark N. Melanson  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K003659