



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 22 2000

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Mr. Michael Turanchik
Director, Research & Development
Medtox Diagnostics, Inc.
1238 Anthony Road
Burlington, North Carolina 27215

Re: K003687
Trade Name: Profile® - II
Regulatory Class: II
Product Code: DKE, DIO, DJG, DKZ, LCM
Dated: November 29, 2000
Received: November 30, 2000

Dear Mr. Turanchik:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.


A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in cursive script that reads "Steven Gutman".

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Attachment 2 (B)

Indications for Use Statement

510(k)
Number
(if known)

K003687

Device Name Verdict®-II THC/COC/OPI/AMP

Indications for Use

The Verdict®-II THC/COC/OPI/AMP test is a one-step immunochromatographic tests for the rapid qualitative detection of: Cannabinoids (THC), cocaine, opiates, and amphetamines in human urine. The cutoffs of these drugs are at the following concentrations:

THC	Tetrahydrocannabinol	50 ng/mL
COC	Cocaine (Benzoylecgonine)	300 ng/mL
OPI	Opiates (Codeine/Morphine)	300 ng/mL
AMP	Amphetamine	1000 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II THC/COC/OPI/AMP test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

Heaven Cooley
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number

K003687

Attachment 2 (C)

Indications for Use Statement

510(k)
Number
(if known)

K 003687

Device Name Verdict® -II THC/COC/OPI

Indications for Use

The Verdict® -II THC/COC/OPI test is a one-step immunochromatographic tests for the rapid qualitative detection of: Cannabinoids (THC), cocaine and opiates in human urine. The cutoffs of these drugs are at the following concentrations:

THC	Tetrahydrocannabinol	50 ng/mL
COC	Cocaine (Benzoylecgonine)	300 ng/mL
OPI	Opiates (Codeine/Morphine)	300 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict® -II THC/COC/OPI test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result


PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐


(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number

K 003687

Attachment 2 (D)

Indications for Use Statement

510(k)
Number
(if known)

K 003687

Device Name Verdict® -II THC/COC/PCP

Indications for Use

The Verdict® -II THC/COC/PCP test is a one-step immunochromatographic tests for the rapid qualitative detection of : Cannabinoids (THC), cocaine and phencyclidines (PCP) in human urine. The cutoffs of these drugs are at the following concentrations:

THC	Tetrahydrocannabinol	50 ng/mL
COC	Cocaine (Benzoylecgonine)	300 ng/mL
PCP	Phencyclidine	25 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict® -II THC/COC/PCP test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801. 109)

OR

Over-The-Counter Use _____

Division Sign-Off)

Division of Clinical Laboratory Devices

Number

K 003687

Attachment 2 (E)

Indications for Use Statement

510(k)

Number

(if known)

K 003687

Device Name Verdict®-II THC/COC/AMP

Indications for Use

The Verdict®-II THC/COC/AMP test is a one-step immunochromatographic tests for the rapid qualitative detection of : Cannabinoids (THC), cocaine and amphetamines in human urine. The cutoffs of these drugs are at the following concentrations:

THC	Tetrahydrocannabinol	50 ng/mL
COC	Cocaine (Benzoylecgonine)	300 ng/mL
AMP	Amphetamine	1000 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II THC/COC/AMP test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒ (Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

Jean Cooper
(Division Sign-Off)

Division of Clinical Laboratory Devices

Number K 003687

Attachment 2 (F)

Indications for Use Statement

510(k)

Number
(if known)

K 003687

Device Name Verdict®-II THC/COC

Indications
for Use

The Verdict®-II THC/COC test is a one-step immunochromatographic tests for the rapid qualitative detection of : Cannabinoids (THC) and cocaine in human urine. The cutoffs of these drugs are at the following concentrations:

THC	Tetrahydrocannabinol	50 ng/mL
COC	Cocaine (Benzoylecgonine)	300 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II THC/COC test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801. 109)

OR

Over-The-Counter Use _____

Jean Cooper
(Division Sign-Off)

Division of Clinical Laboratory Devices

Device Number

K003687

Attachment 2 (G)

Indications for Use Statement

510(k)
Number
(if known)

K 00 3687

Device Name Verdict®-II THC

Indications
for Use

The Verdict®-II THC test is a one-step immunochromatographic tests for the rapid qualitative detection of Cannabinoids (THC in human urine. The cutoff of this drug is at the following concentration:

THC Tetrahydrocannabinol 50 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II THC test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR Over-The-Counter Use _____

Jean Cooper
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K003687

Attachment 2 (H)

Indications for Use Statement

510(k)
Number
(if known)

K 003687

Device Name Verdict®-II COC

Indications for Use

The Verdict®-II COC test is a one-step immunochromatographic tests for the rapid qualitative detection of cocaine in human urine. The cutoff of this drug is at the following concentration:

COC Cocaine (Benzoylecgonine) 300 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II COC test provides only preliminary analytical test results. A more ~~specific alternate chemical~~ method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801. 109)

OR

Over-The-Counter Use _____

Jean Cozger
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K003687

Attachment 2 (I)

Indications for Use Statement

510(k)
Number
(if known)

K003687

Device Name Verdict®-II AMP

Indications for Use

The Verdict®-II AMP test is a one-step immunochromatographic tests for the rapid qualitative detection of amphetamines in human urine. The cutoff of this drug is at the following concentration:

AMP Amphetamine 1000 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II AMP test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

Dean Cooper
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K003687

Attachment 2 (J)

Indications for Use Statement

510(k)
Number
(if known)

K 003687

Device Name Verdict®-II OPI

Indications for Use

The Verdict®-II OPI test is a one-step immunochromatographic tests for the rapid qualitative detection of opiates in human urine. The cutoff of this drug is at the following concentration:

OPI Opiates (Codeine/Morphine) 300 ng/mL

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II OPI test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801. 109)

OR Over-The-Counter Use _____

Jean Cooper
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K 003687

Attachment 2 (K)

Indications for Use Statement

510(k)
Number
(if known)

K 50 3687

Device Name Verdict®-II PCP

Indications for Use

The Verdict®-II PCP test is a one-step immunochromatographic tests for the rapid qualitative detection of phencyclidines (PCP) in human urine. The cutoff of this drug is at the following concentration:

PCP	Phencyclidine	25 ng/mL
-----	---------------	----------

This product is intended for use under medical supervision in hospitals, physician offices, health clinics and drug treatment/counseling centers. It is not for over the counter sale.

The Verdict®-II PCP test provides only preliminary analytical test results. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result


PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801. 109)

OR

Over-The-Counter Use ☐


(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K003687