

K003771

Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim)
Ansell
1875 Harsh Avenue SE
Massillon, Ohio 44646
Telephone: 330-833-2811
Fax: 330-833-6501

JAN 16 2001

Attachment 4

510(k) Summary

Summary

Attached



Signature

James R. Chatterton
Typed Name

December 1, 2000
Date

Ansell

PROFESSIONAL HEALTHCARE

1875 Harsh Ave. S.E. • P.O. Box 550
 Massillon, OH 44648-0550 U.S.A.
 330.833.2811 / 800.321.9752 U.S.A.
 330.833.5991 Fax
 ansellhealthcare.com

Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim)
 Ansell
 1875 Harsh Avenue SE
 Massillon, Ohio 44646
 Telephone: 330-833-2811
 Fax: 330-833-6501

[1] Summary

[2] Ansell
 1875 Harsh Avenue SE
 Massillon, Ohio 44646

Contact: James R. Chatterton
 Telephone: 330-833-2811
 Fax: 330-833-6501

December 1, 2000

- [3] Trade Name: Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim)
 Common Name: Surgical Gloves
 Classification Name: Surgeon's Glove
- [4] Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim) meet all of the requirements of ASTM D 3577-00, Type 1.
- [5] Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim) meet all the current specifications for ASTM D 3577-00 Rubber Surgical Gloves.
- [6] Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim) are sterile disposable devices intended to be worn by operating room personnel to protect a surgical wound from contamination.
- [7] Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim) are summarized with the following technological characteristics compared to ASTM or equivalent standards.

Characteristics	Standard
Dimensions	Meets ASTM D 3577-00
Physical Properties	Meets ASTM D 3577-00, Type 1
Freedom from holes	Meets ASTM D 3577-00 Meets ASTM D 5151-99
Protein Label Claim	This latex glove contains 50 micrograms or less of total water extractable protein per gram.

Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim)
Ansell
1875 Harsh Avenue SE
Massillon, Ohio 44646
Telephone: 330-833-2811
Fax: 330-833-6501

Meets ASTM D 5712-99 Standard Test
Method for Analysis of Protein in Natural
Rubber and Its Products

Biocompatibility

Primary Skin Irritation in Rabbits

Guinea Pig Sensitization

Passes

Passes

- [8] The performance test data of the non clinical tests are the same as mentioned immediately above.
- [9] Clinical data is not needed for medical gloves or for most devices cleared by the 510(k) process.
- [10] It is concluded that the Encore® Orthopaedic Powder Free Latex Surgical Gloves (Protein Label Claim) are as safe, as effective, and perform as well as the glove performance standards referenced above and therefore meet:
- ASTM listed standards,
FDA hole requirements, and
labeling claims for the product.
- [11] This summary will include any other information reasonably deemed necessary by the FDA.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 16 2001

Mr. James R. Chatterton
Vice President Regulatory
Ansell Healthcare Products, Incorporated
1875 Harsh Avenue Southeast
Massillon, Ohio 44646

Re: K003771
Trade Name: Encore Powder-Free Latex Surgical Gloves,
Orthopaedic With Protein Content Labeling Claim
(50 Micrograms or Less)
Regulatory Class: I
Product Code: KGO
Dated: December 19, 2000
Received: December 20, 2000

Dear Mr. Chatterton:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does

not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Attachment 2

Indications for Use Statement

510(k)
Number
(if known)

K003771

Device Name

ENCORE POWDER-FREE LATEX SURGICAL GLOVES, ORTHOPAEDIC WITH PROTEIN
CONTENT LABELING CLAIM (50 MICROGRAMS OR LESS)

Indications for Use

Encore® Orthopaedic Powder Free Latex Surgical Gloves intended use is to be
worn by operating room personnel to protect a surgical wound from contamination.

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED

Concurrence of CDRH Office of Device Evaluation (ODE)

Prescription Use _____
Per 21 CFR 801.109

OR

Over-The-Counter Use _____

David J. Chen

(Division Sign-Off)

Division of Dental, Infection Control,
and General Hospital Devices

510(k) Number K003771