

FEB - 5 2001

K003854

**Infusion Devices, Inc.**  
**730 W. Wilshire, Suite 110**  
**Oklahoma City, OK 73116**

**Non-Confidential Summary of Safety and Effectiveness**

Page 1 of 2

December 12, 2000

Infusion Devices, Inc.  
730 W. Wilshire, Suite 110  
Oklahoma, OK 73116

Tel - (405) 840-4224

Fax - (405) 843-7337

<b>Official Contact:</b>	Clay Kennard - President
<b>Proprietary or Trade Name:</b>	Intravascular Administration Set
<b>Common/Usual Name:</b>	Intravascular Administration Set
<b>Classification Name:</b>	Intravascular Administration Set
<b>Device:</b>	Intravascular Administration Set
<b>Predicate Devices:</b>	Custom Assemblies, Inc. - Intravascular Administration Set – K993463

**Device Description:**

The Infusion Devices Intravascular Administration Set is a tubing set with various connectors, and accessories, such as drip chambers, infusion line filters, clamps, integral check valves, spikes, Y-site injection sites used as an extension of a fluid pathway for I.V. administration.

**Intended Use:**

Indicated Use --	Indicated as a single, use, sterile device for use in I.V. therapy when an extended fluid path is required for administration.
Environment of Use --	Hospital, Emergency Services, Home settings, wherever I.V. administration is utilized.

## Non-Confidential Summary of Safety and Effectiveness

Page 2 of 2

December 12, 2000

### Comparison to Predicate Devices:

Attribute	Infusion Devcies	Custom Assemblies K993463
-----------	---------------------	---------------------------------

### Use

Indicated for I.V. therapy and administration	Yes	Yes
Used as fluid path way	Yes	Yes
Single patient use	Yes	Yes
Environment where I.V therapy is indicated	Yes	Yes

### Design

Tubing of various lengths and diameters	Yes	Yes
Connectors - luer fittings	Yes	Yes
Assembled with - drip chambers, clamps, infusion line filter, spikes, check valves, caps	Yes	Yes
Offered sterile	Yes	Yes

### Materials

Tubing - PVC Class VI	Yes	Yes
Accessory components Polycarbonate, Polypropylene, PVC	Yes	Yes

### Performance Standards / Specifications

None required under Section 514	Yes	Yes
---------------------------------	-----	-----

### Differences between Other Legally Marketed Predicate Devices

There are no significant differences between the intended device and the predicate – Custom Assemblies – K993463.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB - 5 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Infusion Devices, Incorporated  
C/O Mr. Paul E. Dryden  
Consultant  
Promedic, Incorporated  
6329 West Waterview Court  
McCordsville, Indiana 46055

Re: K003854  
Trade Name: Intravascular Administration Set  
Regulatory Class: II  
Product Code: FPA  
Dated: December 12, 2000  
Received: December 13, 2000

Dear Mr. Dryden:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. "A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531


Page 2 - Mr. Dryden

through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K003854

**SECTION 2**

**INDICATIONS FOR USE**

Page 1 of 1

Pursuant to the Notice of February 6, 1996 regarding listing of Indications for Use on a separate sheet, the following is per that request.

**510(k) Number:** \_\_\_\_\_ (to be assigned)

**Device Name:** Intravascular Administration Sets

**Intended Use:** Indicated as a single use, sterile device for use in I.V. therapy when an extended fluid path is required for administration.

**Environment of use:** Hospital, Emergency Medical Services, Home care settings, wherever I.V. fluid administration may be indicated.


---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Prescription Use** ☒   
 (Per CFR 801.109)

or

**Over-the-counter use** ☐

  
(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K003854