
JAN 1 7 2001



X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) Ansell Healthcare Products Inc. 1875 Harsh Avenue SE Massillon, Ohio 44646

Telephone: 330-833-2811

Fax: 330-833-6501

Atta	ch	m	_	nt	A
MILA	GH		u	111	-

		510(k) Summary	
Şummary	Attached		

Signature

James R. Chatterton
Typed Name

December 15, 2000
Date

Ansell
PROFESSIONAL HEALTHCARE

K003895

1875 Harsh Ave. S.E. • P.O. Box 550 Massillon, OH 44648-0550 U.S.A. 330.833.2811 / 800.321.9752 U.S.A. 330.833.5991 Fax ansellhealthcare.com

X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less)

Ansell Healthcare Products Inc.

1875 Harsh Avenue SE

Massillon, Ohio 44646

Telephone: 330-83 Fax: 330-833-6501

[1] Summary

[2] Ansell Healthcare Products Inc. Inc. 1875 Harsh Avenue SE Massillon, Ohio 44646

Contact:

James R. Chatterton

Telephone:

330-833-2811

Fax:

330-833-6501

December 15, 2000

[3] Trade Name:

X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content

Labeling Claim (50 micrograms or less)

330-833-2811

Common Name:

Surgical Gloves

Classification Name:

Surgeon's Glove

- [4] X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) meet all of the requirements of ASTM D 3577-00, Type 1.
- [5] X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) meet all the current specifications for ASTM D 3577-00 Rubber Surgical Gloves.
- [6] X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) are sterile disposable devices intended to be worn by operating room personnel to protect a surgical wound from contamination.
- [7] X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) are summarized with the following technological characteristics compared to ASTM or equivalent standards.

Characteristics

Standard

Dimensions

Meets ASTM D 3577-00

Physical Properties

Meets ASTM D 3577-00, Type 1

Freedom from holes

Meets ASTM D 3577-00

Meets ASTM D 5151-99

K003895

X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) Ansell Healthcare Products Inc. 1875 Harsh Avenue SE Massillon, Ohio 44646

Telephone: 330-833-2811

Fax: 330-833-6501

Powder-Free

Meets ASTM D 3577-00

Not more than 2 mg residue by mass per

glove.

Protein Label Claim

This latex glove contains 50 micrograms or less of total water extractable protein per

gram.

Biocompatability

Primary Skin Irritation in Rabbits

Guinea Pig Sensitization

Passes Passes

- [8] The performance test data of the non clinical tests are the same as mentioned immediately above.
- [9] Clinical data is not needed for medical gloves or for most devices cleared by the 510(k) process.
- [10] It is concluded that the X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) (Modified) are as safe, as effective, and perform as well as the glove performance standards referenced above and therefore meet:

ASTM listed standards, FDA hole requirements, and labeling claims for the product.

[11] This summary will include any other information reasonably deemed necessary by the FDA.



JAN 1 7 2001

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. James R. Chatterton Vice President Regulatory Ansell Healthcare Products, Incorporated 1875 Harsh Avenue, Southeast Massillon, Ohio 44646-7199

Re: K003895

Trade Name: X-Tenda Cuff® Powder Free Latex Surgical

Gloves with Protein Content Labeling Claim

(50 Micrograms or Less)

Regulatory Class: I Product Code: KGO

Dated: December 15, 2000 Received: December 18, 2000

Dear Mr. Chatterton:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does

not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address

"http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours

Timothy A. Ulatowski

Director

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less)

Ansell Healthcare Products Inc.

1875 Harsh Avenue SE Massillon, Ohio 44646

Telephone:

330-833-2811

Fax:

330-833-6501

Attachment 2

Indications for Use Statement

510(k) Number (if known)	K003895
Device Name	X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less)
Indications for Use	X-Tenda Cuff® Powder Free Latex Surgical Gloves with Protein Content Labeling Claim (50 micrograms or less) intended use is to be worn by operating

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH Office of Device Evaluation (ODE)

/Division Sign-Off)
vision of Dental, Infection Control,

Prescription Use Per 21 CFR 801.109

Over-The-Counter Use X