

FEB - 9 2001

K010095

**APPENDIX I**

**510(k) SUMMARY**

SUMMARY OF THE SAFETY AND EFFECTIVENESS FOR NON-STERILE PURPLE POWDERED LATEX EXAMINATION GLOVES WITH AND WITHOUT GRAPE SCENTING AND WITH A PROTEIN LABELING CLAIM

**Contact person :** Cheah Chor Hee

This summary of safety and effectiveness information is being submitted in accordance with the requirement of SMDA 1990.

**Device Information:**

Trade Name – NON-STERILE PURPLE POWDERED LATEX EXAMINATION GLOVES WITH & WITHOUT GRAPE SCENTING.

Common Name - Exam gloves

Classification Name - Patient examination glove (per 21 CFR 880.6250)

Classification Information - Class I latex patient examination glove 80LYY, powdered and meeting all the requirements of ASTM-D3578-99 Standard Specification for Latex Examination Gloves for Medical Application.

**Device Description:**

Class I latex patient examination gloves 80LYY, powdered and meeting all the requirements of ASTM-D3578-99 Standard Specification for Latex Examination Gloves for Medical Application.

**Intended Use of Device:**

A medical glove to be worn on the hand of the health care and similar personnel to prevent contamination between health care personnel and patient.

Technological Characteristics of Device:

**1. Dimension**

DIMENSION					
		Ambidextrous		Size Fitted	
Width		X-Small	70 mm +/- 10 mm	5.5	70 +/- 10 mm
		Small	80 mm +/- 10mm	6.0	76 +/- 10mm
		Medium	95 mm +/- 10mm	6.5	83 +/- 10mm
		Large	111mm +/- 10mm	7.0	89 +/- 10mm
				7.5	95 +/- 10mm
				8.0	102 +/- 10mm
				8.5	108 +/- 10 mm
				9.0	114 +/- 10mm
Length		230 mm min			
Thickness -	Finger	0.08 mm min			
	Palm	0.08 mm min			

## 2. Physical Properties (ASTM-D3578-99 Standard Specification for Latex Exam Gloves)

LOT #	TENSILE STRENGTH				ULTIMATE ELONGATION			
	AGED		UNAGED		AGED		UNAGED	
TESTED	SGMP	ASTM	SGMP	ASTM	SGMP	ASTM	SGMP	ASTM
X-SMALL 0178	27.2	14.0	27.3	14.0	870	500	900	700
SMALL 0178	28.8	14.0	28.0	14.0	900	500	930	700
MEDIUM 0178	27.5	14.0	26.1	14.0	870	500	910	700
LARGE 0178	26.6	14.0	27.8	14.0	840	500	910	700

## 3. Water Tight Test Data

BATCH NUMBER	DATE TESTED	SAMPLING SIZE	LEAK STATUS	NUMBER LEAKED
<b>Unaged Smpl</b>				
0178 XS	6 Dec 00	125	Yes	1
0178 S	6 Dec 00	125	Yes	2
0178 M	6 Dec 00	125	No	0
0178 L	6 Dec 00	125	Yes	1
<b>Aged Smpl</b>				
0178 XS	15 Dec 00	125	Yes	1
0178 S	15 Dec 00	125	Yes	1
0178 M	15 Dec 00	125	No	0
0178 L	15 Dec 00	125	Yes	2

The above figures are within the ASTM D-3578-99 requirements for latex exam gloves of 2.5% AQL.

## 4. Biocompatibility

### BIOCOMPATIBILITY TESTS

Test results indicate that the gloves passed the biocompatibility tests for gloves.

## 5. Residual Protein Level

TESTS	FDA ALLOWABLE LEVEL	CLAIMED LEVEL
ASTM D 5712-95	-	< 200 µg/g Range: 76 - 118 µg/g Mean: 95 µg/g

The data presented indicates that the Non-sterile Latex Powdered latex examination glove

1. meets/exceeds ASTM- D3578-99 Standard Specifications For Latex Examination Glove,
2. meets FDA pinhole requirements,
3. meets the protein labeling claim level at  $<200 \mu\text{g/g}$ .



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB - 9 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SGMP Company Limited  
C/O Ms. Janna P. Tucker  
Official Correspondent for  
Tucker & Associates  
198 Avenue De La D' Emerald  
Sparks, Nevada 89434

Re: K010095  
Trade Name: Nonsterile Powdered Purple Latex  
Examination Gloves, with & without Grape Scent with  
Protein Labeling (200 Micrograms or Less)  
Regulatory Class: I  
Product Code: LYY  
Dated: January 8, 2001  
Received: January 11, 2001

Dear Ms. Tucker:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



*ty* Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**INDICATION FOR USE STATEMENT**

Applicant : SGMP Company Limited

510K NUMBER : K010095

Device Name : Powdered Latex Examination Gloves, PURPLE WITH/WITHOUT GRAPE SCENT  
AND WITH A PROTEIN LABELING CLAIM (ZOOMERLESS)

Indication For Use :

This is a medical glove to be worn on the hand of health care and similar personnel to prevent contamination between health care personnel and the patient.

.....  
Concurrence of CDRH Office of Device Evaluation (ODE)

Prescription Use .....  
Per 21 CFR 801.109

OR

Over-The-Counter...X.....

2

Chia S. Lim  
(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices  
510K Number K010095