

FEB 14 2001

K010212

SUMMARY OF SAFETY AND EFFECTIVENESS

Sponsor: Biomet Inc.
Airport Industrial Park
P.O. Box 587
Warsaw, IN 46581-0587

Contact Person: Carol Lauster

Device(s): stem adapter assembly, adapter, locking insert

Classification: Class II

Device Product Code: 87 JDW (21 CFR 888.3690)

Intended Use: The offset tibial tray is for cemented use and has the same indications as MCK knee system, 1) Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved, 2) Correction of valgus, varus or posttraumatic deformity, 3) Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure.

These devices are for cemented use only.

Device Description: The offset tibial tray is designed with sizes 0mm, 2.5mm and 5mm, stem adapters. The offset tray system will be used with the existing stem extensions and is connected with a morse taper as is used in standard revision trays. The stem, adapter, and tray will lock together via taper and screw mechanism. The stem adapter assembly, adapter and locking insert will be packaged together.

The offset tibial tray may be used with the Maxim (K915132), Ascent (K982869) and AGC (K833921) Knee Systems.

Materials: Offset tibial tray components are made from Ti-6Al-4V (ASTM F136)

Potential Risks: The potential risks associated with this device are the same as with any joint replacement device. These include, but are not limited to:

Reaction to bone cement	Blood vessel damage	Bone fracture
Deformity of the joint	Soft tissue imbalance	Infection
Cardiovascular disease	Delayed wound healing	Hematoma
Fracture of the cement	Metal sensitivity	Dislocation
Implant loosening/migration	Fracture of the components	Excessive wear
Tissue growth failure	Nerve damage	Valgus-varus deformity

Substantial Equivalence: MCK Knee System, K915132

000020



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB 14 2001

Ms. Carol Lauster
Regulatory Specialist
Biomet, Inc.
P.O. Box 587
Warsaw, Indiana 46581

Re: K010212
Trade Name: Offset Tibial Tray
Regulatory Class: II
Product Code: MBV
Dated: January 22, 2001
Received: January 23, 2001

Dear Ms. Lauster:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Carol Lauster

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a long horizontal flourish extending to the right.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K010212
Device Name: Offset Tibial Tray

Indications for Use: 1) Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved. 2) Correction of varus, valgus, or posttraumatic deformity. 3) Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure

These devices are for cemented use only.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____
(Optional Format 1-2-96)

for Mark A. Melkerson
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K010212