

K010284 1/2

MAY - 1 2001

**510(k) Summary Statement  
For the Lyra G™ Surgical Laser System, Accessories &  
Angled Delivery Devices for the Treatment of Benign Prostatic Hyperplasia (BPH)  
and Pseudofolliculitis barbae (PFB)**

---

**General Information**

A. Trade Name  
Lyra G™ Series Surgical Laser System (SL Series Q-Switched Nd:YAG configuration)

B. Common Name  
Laser Instrument, Surgical, Powered

C. Establishment Registration Number  
  
2937094

D. Manufacturer's Identification  
  
Laserscope  
3070 Orchard Drive  
San Jose, CA 95134-2011  
(408) 943-0636  
(503) 961-1688 FAX  
  
Official Correspondent  
Paul Hardiman  
Manager, Regulatory Affairs/Clinical Affairs

E. Device Classification  
  
The Lyra G Series Surgical Laser System has been specifically classified as a Class II medical device by the OB/GYN, General Plastic Surgery, and ENT Device Advisory Panels.

F. Performance Standards  
  
The Lyra G Series Surgical Laser System conforms with federal regulations and the performance standards 21 CFR 1040.10 and 1040.11 for medical laser systems.

G. Predicate Devices:

- Laserscope Lyra™ Laser System and Accessories
- Laserscope Lyra G™ Laser System and Accessories
- Modified Coherent VersaPulse Select Single Wavelength (Ho:YAG) and Dual Wavelength (Ho:YAG/Nd:YAG) Surgical Lasers and Delivery Devices with Accessories

**H. Product Description:**

The Laserscope Lyra G Surgical Laser System, Accessories, Angled Delivery Devices and Cooling Device are comprised of the following main components:

- A Laser Console
- A Fiber Port (for Delivery Devices)
- Control and Display Panels
- Operating Software
- Footswitch and Handswitch Delivery Controls
- A variety of Delivery Devices and Accessories
- A Cooling Sub-system

**I. Indications For Use:**

The Lyra G™ Series Surgical Laser System, Accessories, and Angled Delivery Devices are intended for use in cutting, coagulating, and vaporizing of prostatic tissue during treatment of Benign Prostatic Hyperplasia ( BPH) with the KTP wavelength @ 532nm. The device is not intended to treat prostate cancer. When used at 1064nm the Lyra G™ is intended to be used for the removal of unwanted body hair which causes the condition known as Pseudofolliculitis barbae (PFB) (Shaving/Razor Bumps) in Fitzpatrick Skin Types I to VI.

**J. Rationale for Substantial Equivalence**

The Laserscope the Lyra G™ Surgical Laser System, Accessories/Delivery devices, and Cooling Device share the same indications for use, similar design features, functional features, and therefore are substantially equivalent to: Laserscope Lyra Laser System and Accessories; and, the Modified Coherent VersaPulse Select Wavelength (Ho:YAG) and Dual Wavelength (Ho:YAG/Nd:YAG) Surgical Lasers and Delivery Devices and Accessories. Details are provided in the Substantial Equivalence Section of this submission.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY - 1 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Paul H. Hardiman  
Manager, Regulatory Affairs/Clinical Affairs  
Laserscope  
3070 Orchard Drive  
San Jose, California 95134

Re: K010284  
Trade/Device Name: Lyra G™ Series Surgical Laser System  
Regulation Number: 878.4810  
Regulatory Class: II  
Product Code: GEX  
Dated: January 25, 2001  
Received: January 31, 2001

Dear Mr. Hardiman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

*Miriam C. Probst*  
for Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

Page 1

510(k) Number:

K010284

Device Name:

LYRA G™ SERIES SURGICAL LASER SYSTEM

### INTENDED USE:

The Laserscope Lyra G™ laser system has been developed for selective cutting and coagulating of soft tissue. More specifically when used at 532 nm it is intended to hemostatically ~~vaporize prostate tissue~~ of men suffering from benign prostate hyperplasia/hypoplasia (BPH). The device is not intended to treat prostate cancer. When used at 1064 nm the Lyra G is intended to be used for the removal of unwanted body hair which causes the condition known as Pseudofolliculitis barbae (PFB) (Shaving/Razor Bumps) in Fitzpatrick Skin Types I to VI.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: ✓

or

Over The-Counter-Use

(per 21 CFR 801.109)

Miriam C. Probst

(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

510(k) Number K010284

000005