

MAY 10 2001

BioHorizons Implant Systems, Inc.
510(K) Notification
February 15, 2001

K010458

510(k) Summary of Safety and Effectiveness

Proprietary Name

The Maestro System™

Common Name

Uncoated and HA-coated screw-form implants

Classification Name

Endosseous implants, surgical components, and prosthetic attachments

Classification

Class III

Official Contact

*R. Steven Boggan, M.S., M.B.A.
President and Chief Executive Officer
BioHorizons Implant Systems, Inc.
One Perimeter Park South
Suite 230 South
Birmingham, AL 35243
(205) 967-7880
FAX (205) 870-0304*

Device Description

The Maestro System™ is a comprehensive system containing implants, surgical components, and prosthetic components. The implants are specifically designed to optimize strain distribution to contiguous bone under functional loading in order to promote strain-induced bone growth and interface maintenance over the long-term. This improvement in biomechanical performance is achieved by optimizing implant designs specifically for each bone density classification (D1, D2, D3, and D4) and bone volume classification (Division A, B, and C-h) in the mandible and maxilla.

Three base implant designs, corresponding to each bone density classification (D1/D2, D3 and D4), are available in $\phi 3.5$ mm, $\phi 4.0$ mm and $\phi 5.0$ mm diameters. Each implant design, manufactured from titanium alloy conforming to ASTM F 136, is available in two lengths and may feature a resorbable blast media (RBM) treated surface, or hydroxylapatite (HA) coating. The following table provides a comprehensive summary of the proposed implant sizes.

Diameter	Design	Lengths (mm)	Coating
$\phi 3.5$	D3	9, 12	RBM
$\phi 4.0$	D2	9, 12	RBM
	D3	9, 12	RBM
	D4	9, 12	HA
$\phi 5.0$	D2	9, 12	RBM
	D3	9, 12	RBM
	D4	9, 12	HA

Table. Proposed implant sizes.

Product Evaluation

Evaluation of The Maestro System™ consisted of mechanical testing of the implant and bioactive coating mechanical tests. These analyses indicate The Maestro System™ should be safe and effective when used as intended.

Indications

The Maestro System™ may be used in the mandible and maxilla for use as an artificial root structure for single tooth replacement or as abutments for fixed bridgework and denture retention.

Substantial Equivalence Information

The Maestro System™ is substantially equivalent in all features which could affect safety or effectiveness to the BioHorizons Dental Implant System (K96433, K964330 and K972313).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY 10 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Steven Boggan
President of Operations
Biohorizons Implant Systems, Incorporated
One Perimeter Park South, Suite 230 South
Birmingham, Alabama 35243

Re: K010458
Trade/Device Name: The Maestro System
Regulation Number: 872.3640
Regulatory Class: III
Product Code: DZE
Dated: February 15, 2001
Received: February 16, 2001

Dear Mr. Boggan:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note:

this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K010458


Device Name: The Maestro System

Indications for Use:

The Maestro SystemTM may be used in the mandible and maxilla for use as an artificial root structure for single tooth replacement or as abutments for fixed bridgework and denture retention.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K010458

Prescription Use ☒ ✓
(per 21 CFR 801.109)

OR

Over-The-Counter Use