

AUG 1 0 2001

510K #K010544 Submission for Vertex L/C Orthodontic Direct Bonding Paste SH 12 of 15
Apex Dental Materials, 603 Berkley Court
Schaumburg, IL. 60194

510 (K) SUMMARY

As Required by the Safe Medical Devices Act of 1990

Apex Dental Materials, Inc.
603 Berkley Court
Schaumburg, IL. 60194
Phone: (847) 490-1014

510 (K) Submission Date: February 22, 2001

Contact Person: Scott Lamerand

Device Name:

Trade Name:	Vertex L/C Orthodontic Direct Bonding Paste
Common Name:	Orthodontic Bracket Adhesive
Classification Name:	Bracket Adhesive Resin and Tooth Conditioner, per 21 CFR parts 872.3750

Classification:

Regulatory Class:	II
Product Code:	DYH

IDENTIFICATION OF THE LEGALLY MARKETED PREDICATE DEVICE

PREDICATE DEVICE

Light Bond Bonding Resin™

Light Bond Bonding Paste™ (Reliance Orthodontic Products) is a single paste composite orthodontic bonding system, utilizing a BisGMA based resin paste between a properly prepared etched enamel tooth surface with sealant, and the orthodontic bracket. Its physical properties are similar to the applicant device and uses are identical. Like the applicant device, Light Bond Bonding Paste™ is one component used in conjunction with a complete orthodontic bonding system. It hardens by a light cure polymerization mechanism employing a light initiator, and a chemical activator.

Summary continued:

DESCRIPTION OF APPLICATION DEVICE

VERTEX L/C ORTHODONTIC DIRECT BONDING PASTE

Vertex L/C Orthodontic Direct Bonding Paste is designed to provide an orthodontist with a high strength, rapid cure, resin based paste material for use in improving the retention of metal, ceramic (and/or plastic) orthodontic brackets to etched-enamel tooth surfaces. The product is a BisGMA based resin paste used between a properly prepared orthodontic sealant and a metal, ceramic (and/or plastic) orthodontic bracket. The product is cured via photo-initiated free radical polymerization. When properly employed and used with an orthodontic sealant to an etched enamel tooth surface, the paste is designed to maintain bracket adhesion for the duration of the orthodontic treatment. Vertex L/C Orthodontic Direct Bonding Paste is designed to be marketed as a stand alone product but can be made available in kit form including Vertex Etchant (K010849) and Vertex Sealant L/C Orthodontic Sealant (K010547).

INTENDED USES OF APPLICANT DEVICE

Vertex L/C Orthodontic Direct Bonding Paste is indicated to provide adhesion between an etched enamel substrate tooth surface and a metal, ceramic (and/ or plastic) orthodontic bracket when using a bonding sealant, such as Vertex Sealant L/C Orthodontic Sealant (K010547).

Summary continued:

PERFORMANCE CHARACTERISTICS and CONCEPTS

Vertex L/C Orthodontic Direct Bonding Paste has similar handling to the Light Bond Bonding Paste™. From the physical testing observations and analysis, including shear bond strength, diametral tensile strength and compressive strength, we suggest that Vertex L/C Paste is substantially equivalent to Light Bond Bonding Paste™ (Reliance Orthodontic Products). Along with this we would suggest the individual components of Vertex L/C Paste are long time industry standards and are utilized in numerous orthodontic bracket-bonding systems currently marketed in the United States (see Confidential Formulation Details on page 5).

Equivalent Product and Manufacturer

Corresponding 510(k) Numbers

Light Bond™ (Reliance Orthodontic Products)

K880793



AUG 1 0 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Scott Lamerand
Apex Dental Materials, Incorporated
603 Berkley Court
Schumburg, Illinois 60194

Re: K010544
Trade/Device Name: Vertex L/C Orthodontic Director
Bonding Paste
Regulation Number: 872.3750
Regulatory Class: II
Product Code: DYH
Dated: August 3, 2001
Received: August 3, 2001

Dear Mr. Lamerand:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director

Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K010544

510K #K010544 Submission for Vertex L/C Orthodontic Direct Bonding Paste SH 4 of 15
Apex Dental Materials, 603 Berkley Court
Schaumburg, IL. 60194

Indications for Use

510(K) Number (if known): K010544

Device name: Vertex L/C Orthodontic Direct Bonding Paste

Indications For Use:

Vertex L/C Orthodontic Direct Bonding Paste is designed to provide an orthodontist with a high strength, rapid cure, resin based paste material for use in improving the retention of metal, ceramic (and/or plastic) orthodontic brackets to etched-enamel tooth surfaces. The product is a BisGMA based resin paste used between a properly prepared orthodontic sealant and a metal, ceramic (and/or plastic) orthodontic bracket. The product is cured via photo-initiated free radical polymerization. When properly employed and used with an orthodontic sealant to an etched enamel tooth surface, the paste is designed to maintain bracket adhesion for the duration of the orthodontic treatment.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)


Concurrence of CDRL, Office of Device Evaluation (ODE)

Prescription use X
(Per 21 CFR 801.109)

OR

Over- The- Counter Use

(Optional Format 1-2-96)


(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K010544