

1K010548

MAR 12 2001

510(k) Summary

This summary of the 510(k) safety and effectiveness information is submitted in accordance with the requirements of the Safe Medical Devices Act of 1990 and 21 C.F.R. 807.92.

- 1 The submitter of this premarket notification is:

Kevin J. O'Connell
Regulatory Engineer
Radionics, a division of Tyco Healthcare Group LP
22 Terry Avenue
Burlington, MA 01803
Tel.: (781) 272-1233
Fax: (781) 272-2428

This summary was prepared on Tuesday, February 20, 2001.

- 2 The name of the device is the Radionics Screw-In Dynamic Reference Frame (SI-DRF) for use with Radionics Optical Tracking System (OTS). The common name is an Intraoperative Guidance Device, and its classification name is a stereotaxic instrument (accessory).
- 3 The above device is substantially equivalent to the Radionics Spinal DRF for OTS and Medtronic StealthStation's Fixation Tube.
- 4 The above device consists of a light array and shaft that attaches to a commercially available distraction post. When coupled with the OTS workstation, the device allows for preoperative and operative planning of surgical procedures through workstation images.
- 5 The device is similar to its predicates, in that they are all used to attach an array to the patient to maintain registration of the system during surgical procedures.
- 6 The technological characteristics are the same or similar to those found with the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAR 12 2001

Mr. Kevin J. O'Connell
Senior Regulatory Associate
Radionics, a division of Tyco Healthcare Group LP
22 Terry Avenue
Burlington, Massachusetts 01803

Re: K010548
Trade Name: Radionics Screw-In Dynamic Reference Frame
Regulatory Class: II
Product Code: HAW
Dated: February 23, 2001
Received: February 26, 2001

Dear Mr. O'Connell:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Kevin J. O'Connell

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for Celia M. Witten, Ph.D., M.D.
Director

Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K010548

Device Name: Optical Tracking System (OTS) with Screw-In DRF

Indications for use:

The OTS is a graphical planning tool that allows for preoperative planning of, and intraoperative navigation during, cranial and spinal surgical procedures. This includes use in frameless cranial biopsy procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

PRESCRIPTION USE ☒

OR

Over-The-Counter Use ☐

(Per 21 CFR 801.109)

(Optional Format 1-2-96)

for Mark N. Milken
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K010548