

SEP 1 2 2001

K010602

510 (k) Summary of Safety and Effectiveness for VectorVision Hip

Manufacturer:

Address: BrainLAB AG
Ammerthalstrasse 8
85551 Heimstetten
Germany
Phone: +49 89 99 15 68 0
Fax: +49 89 99 15 68 33

Contact Person: Mr. Rainer Birkenbach

Summary Date: February 23, 2001

Device Name:

Trade name: VectorVision Hip

Common/Classification Name: VectorVision, BrainLAB Image Guided Surgery System / Instrument,
Stereotaxic

Predicate Device:

Vector Vision² (K 983831)

Device Classification Name: Instrument, Stereotaxic

Regulatory Class: Class II

Intended Use:

BrainLAB VectorVision is intended to be an intraoperative image guided localization system to enable minimally invasive surgery. It links a freehand probe, tracked by a passive marker sensor system to virtual computer image space on a patient's preoperative image data being processed by a VectorVision workstation. The system is indicated for any medical condition in which the use of stereotactic surgery may be appropriate and where a reference to a rigid anatomical structure, such as the skull, a long bone, or vertebra, can be identified relative to a CT, X-ray or MR based model of the anatomy.

Example procedures include but are not limited to:

Orthopedic and Reconstructive Surgery such as:

Total Joint Replacement (TJR)

Revision surgery of TJR

Minimal Invasive Orthopedic Surgery

Tumor Resection and bone/joint reconstruction

Device Description:

BrainLAB VectorVision Hip is intended to be an intra-operative image guided localization system to perform orthopedic surgery. It contains a pre-operative spatial planning in particular for total joint replacement according to semi-automatically detected anatomical landmarks and axes. During the planning process commercially available product lines of implant manufacturers can be displayed in the patient's anatomy and fine tuned with regard to their seizing and positioning. VectorVision Hip allows semi-automatic correlation of preoperative image data and planning to the actual patient. The software links a freehand probe, tracked by a passive marker sensor system to virtual computer image space on a patient's preoperative image data being processed by a VectorVision workstation. Furthermore it allows calibration and integration of surgical tools to the navigation process.

Substantial equivalence:

VectorVision Hip has been verified and validated according to BrainLAB's procedures for product design and development. The validation proves the safety and effectiveness of the system. The information provided by BrainLAB in this 510 (k) application was found to be substantially equivalent with the predicate devices VectorVision (K983831).



SEP 12 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Rainer Birkenbach
Executive Vice President
BrainLAB, AG
Ammerthalstrasse 8
85551 Heimstetten
Germany

Re: K010602
Trade/Device Name: VectorVision Hip
Regulation Number: 882.4560
Regulatory Class: II
Product Code: HAW
Dated: August 21, 2001
Received: August 24, 2001

Dear Mr. Birkenbach:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health



Enclosure

510(k) Number (if known): K010602Device Name: VectorVision Hip**Indications For Use:**

BrainLAB VectorVision is intended to be an intraoperative image guided localization system to enable minimally invasive surgery. It links a freehand probe, tracked by a passive marker sensor system to virtual computer image space on a patient's preoperative image data being processed by a VectorVision workstation. The system is indicated for any medical condition in which the use of stereotactic surgery may be appropriate and where a reference to a rigid anatomical structure, such as the skull, a long bone, or vertebra, can be identified relative to a CT, X-ray or MR based model of the anatomy.

Example procedures include but are not limited to:

Orthopedic and Reconstructive Surgery such as:

Total Joint Replacement (TJR)

Revision surgery of TJR

Minimal Invasive Orthopedic Surgery

Tumor Resection and bone/joint reconstruction

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)


Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format I-2-96)


(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K010602