

MAY 1 8 2001

K010644

29 Wells Avenue, Yonkers, New York 10701

(914) 963-2040

(800) 431-2972

FAX (914) 963-2567

**Tab F: 510(k) Summary of Safety and Effectiveness**Name, address, phone and fax numbers for person submitting the 510(k) notification:

Arnold Silverman, President  
Skil-Care Corporation  
29 Wells Avenue  
Yonkers, NY 10701  
Phone: 914-963-2040  
Fax: 914-963-2567

Contact person: Arnold SilvermanDate summary was prepared: February 14, 2001Device name:

Trade name: Cinch-Type Disposable Limb Holder  
Common name: Limb Holder  
Classification name: Protective Restraint

Predicate device:

*Economy Limb Holder*, Model # 2510, marketed by the J. T. Posey Company.

Device description:

The Cinch-Type Disposable Limb Holder has a 3-inch-wide by 9 1/2-inch long cuff made from Vel-Foam, a 1/4-inch thick polyurethane foam material backed with nylon "loop" Velcro fabric. A one-inch-wide polyester strap 52 inches in length is stitched to the fabric side of the cuff in two places: At the end of the cuff, near the rounded edge and at the center of the cuff. A single Delrin d-ring, 1 1/4 inches wide, is stitched to the cuff at the same point as the webbing positioned at the center of the cuff. A Velcro "hook" tab measuring 1 1/4" x 1 1/4" is stitched to the foam side of the rounded edge of the cuff --the side opposite the strap. The "hook" material secures to the "loop" material on the cuff thereby permitting a caregiver to position the cuff. The polyester strap is then wrapped around the wrist or ankle, threaded through the d-ring (where it is tied with a single loop-through knot) and then secured to the bed frame. The device is sold in pairs.

Summary, Page 2.

Indications for use:

Cinch-Type Disposable Limb Holders are designed for patients who

- \* Disrupt medical treatments such as IVs or catheters
- \* Are prone to self-injury

Comparative information:

The device submitted herein is substantially equivalent to the *Economy Limb Holder*, Model # 2510 legally marketed by the J.T. Posey Company

Note:

The use of all patient restraints in nursing home are subject to Health Care Financing Administration's regulations which prohibit the use of any physical restraint imposed for the purpose of discipline or convenience. Further, most health care facilities are accredited. HCFA rules governing appropriate use and accreditation standards for device use and personnel training provide the control necessary to ensure that the devices are used correctly. The application of these standards along with public awareness and health care provider training have contributed significantly to ensuring that the least restrictive restraint is used, that restraints are used only when needed for proper medical treatment, and that their use is under appropriate supervision.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY 18 2001

Mr. Arnold Silverman  
President  
Skil-Care Corporation  
29 Wells Avenue  
Yonkers, New York 10701

Re: K010644  
Trade/Device Name: Cinch-Type Disposable Limb Holders,  
Model 912025  
Regulation Number: 880.6760  
Regulatory Class: I  
Product Code: FMQ  
Dated: March 5, 2001  
Received: March 5, 2001

Dear Mr. Silverman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

  
for Timothy A. Ulatowski

Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known) Not known K010644

Device Name: Cinch-Type Disposable Limb Holder

Indications for Use:

The Cinch-Type Disposable Limb Holder is intended for patients who

- \* Disrupt medical treatment such as IVs or catheters
- \* Are prone to self-injury

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE  
IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓

-OR-

Over-The Counter Use \_\_\_\_\_

*Patricia Cuente*

(Division Sign-Off)

Division of Dental, Infection Control,

and General Hospital Devices

510(k) Number K010644 -14-

(Optional Format 1-2-96)