

MAY 3 0 2001

**510(k) Summary**

Device Proprietary Name: OsteoMed 2.0/2.4 Cannulated Screw System

Device Common Name: Bone Screw

Classification Name: HWC, Screw, Fixation, Bone

Name of Submitter: OsteoMed Corporation
3750 Realty Road
Addison, Texas 75001
Phone: (972) 241-3401
Fax: (972) 241-3449

Contact Person: Dawn T. Holdeman

Date Prepared: March 9, 2001

Summary:

This submission describes the OsteoMed 2.0/2.4 Cannulated Screw System indicated for bone fixation following trauma or osteotomy. Screws are intended for single use only.

The OsteoMed 2.0/2.4 Cannulated Screw System is comprised of screws in diameters of 2.0mm to 2.4mm in lengths of 6.0mm to 30.0mm. The screws are made from titanium alloy. Depth gauges, screwdrivers, countersinks and preparation instruments will also be a part of the system.

Equivalence for this device is based on similarities in intended use, material, design and operational principle to the OsteoMed Super Screw Fixation System and the Vilex Cannulated Bone Screw.

Due to the similarity of materials and design to both pre-enactment and post-enactment devices, OsteoMed believes that the OsteoMed 2.0/2.4 Cannulated Screw System does not raise any new safety or effectiveness issues.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY 30 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Dawn T. Holdeman
Regulatory Affairs and Document Control
OsteoMed Corporation
3750 Realty Road
Addison, Texas 75001

Re: K010783

Trade/Device Name: OsteoMed 2.0/2.4 Cannulated Screw System
Regulation Number: 888.3040
Regulatory Class: II
Product Code: HWC
Dated: March 14, 2001
Received: March 15, 2001

Dear Ms. Holdeman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Dawn T. Holdeman

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten".

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health


Enclosure

OsteoMed "Indications for Use" Submission

510(k) Number: K010783

Device Name:	OsteoMed 2.0/2.4 Cannulated Screw System
Indications for Use:	<p>Indicated for bone fixation of hand and foot following trauma or osteotomy. Screws are intended for single use only.</p> <p>The OsteoMed 2.0/2.4 Cannulated Screw System is not intended for use in and is contraindicated for: in cases of active or suspected infection or in patients who are immunocompromised; in patients previously sensitized to titanium; or in patients with certain metabolic diseases; in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation.</p>

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K010783

Prescription Use X
(Per 21 CFR 810.109)

Over-The Counter-Use _____
(Optical Format 1-)