

APR - 6 2001

K010801

**Special 510(k) Summary of Safety and Effectiveness: Device Modifications to the Osteo IC Retrograde/Antegrade Nail****Submission Information**

Name and Address of the Sponsor of the 510(k) Submission: Howmedica Osteonics Corp  
59 Route 17  
Allendale, NJ 07401-1677

Contact Person: Karen Ariemma  
Regulatory Affairs Specialist

Date of Summary Preparation: March 15, 2001

**Device Identification**

Proprietary Name: T2 Femoral Nail (formerly the Osteo IC R/A Nail)  
Common Name: Intramedullary Nail, Femoral Nail  
Classification Name and Reference: Intramedullary Fixation Rod, 21 CFR §888.3020

**Predicate Device Identification**

The Osteo IC Retrograde/Antegrade (R/A) Femoral Nail (to be renamed the T2 Femoral Nail System) is a cylindrical, cannulated titanium alloy tube, slightly bowed to accommodate the shape of the femur. The Osteo IC R/A Femoral Nail may be inserted into the femoral canal using either a retrograde or antegrade surgical approach.

**Description of Device Modification**

The design change involves changing the inner hole and the outer diameter of the connection part of the nails Ø9-Ø11mm, the thread length and slot pattern of the Osteo IC A/R Nail to allow for the Advanced Locking Mode used in the T2 Tibial Nails. The Advanced Locking Mode needs a new short compression screw similar to the T2 Tibial Nail System. The end caps, and Locking Screws of the T2 Tibial Nail System will be compatible with the T2 Femoral Nail System.

**Intended Use**

The subject T2 Femoral Nail System, like the predicate Osteo IC R/A Femoral Nail System, is a fracture fixation device comprised of femoral nails and the related locking screws, compression screws, and end caps. The subject and predicate devices are intended to provide strong and stable internal fracture fixation with minimal soft tissue irritation. This device is utilized as an aid to healing, not as a substitute for normal intact bone and tissue.

**Statement of Technological Comparison:**

FEA analysis demonstrates the comparable mechanical properties of the subject T2 Femoral Nail System to the predicate Osteo IC R/A Femoral Nail System.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

APR - 6 2001

Ms. Karen Ariemma  
Regulatory Affairs Specialist  
Howmedica Osteonics Corporation  
59 Route 17  
Allendale, New Jersey 07401-1677

Re: K010801  
Trade Name: T2 Femoral Nail  
Regulatory Class: II  
Product Codes: HSB  
Dated: March 15, 2001  
Received: March 16, 2001

Dear Ms. Ariemma:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 – Ms. Karen Ariemma

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

for *Celia M. Witten*

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and  
Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K 010801

Device Name: T2 Femoral Nail System

Indications for Use

The T2 Femoral Nail is indicated for long bone fracture fixation specifically femoral fracture fixation, which may include the following:

- Open and closed femoral fractures
- Pseudoarthrosis and correction osteotomy
- Pathologic fractures, impending pathologic fractures, and tumor resections
- Supracondylar fractures, including those with intra-articular extension
- Ipsilateral femur fractures
- Fractures proximal to a total knee arthroplasty
- Fractures distal to a hip joint
- Nonunions and malunions

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒

OR

Over-The-Counter Use ☐

(Per 21 CFR 801.109)

(Optional Format 1-2-96)

for Mark A. Millerson  
(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

510(k) Number K 010801