K010862

## XIV. SUMMARY OF SAFETY AND EFFECTIVENESS



## 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS ULTRAFREE MAX STERILE LATEX POWDER-FREE SURGICAL GLOVES

Manufacturer:

Allegiance Healthcare Sdn. Bhd.

Plot 87, Kampung Jawa 11900 Bayan Lepas Penang, West Malaysia

Regulatory Affairs Contact: Erica Sethi

Allegiance Healthcare Corporation 1500 Waukegan Road, MP-WM

McGaw Park, IL 60085

Telephone:

(847) 785-3337

Date Summary Prepared:

3/5/01

Product Trade Name:

Ultrafree Max Sterile Latex Powder-Free Surgical Gloves

Common Name:

Surgical Glove

Classification:

Glove, Surgical

Predicate Devices:

Ultrafree Max Sterile Latex Powder-Free Surgical Gloves

Description:

Ultrafree Max Powder-Free Surgical gloves are formulated

using natural rubber latex and offered sterile.

Intended Use:

Ultrafree Max Sterile Latex Powder-Free Surgical Gloves are intended to be worn by operating room personnel to protect a surgical wound from contamination. These latex gloves contain 50 micrograms or less of total water extractable protein per gram. In addition, these gloves have been tested for use with chemotherapy drugs.

Substantial Equivalence:

Ultrafree Max Sterile Latex Powder-Free Surgical Gloves with protein and chemotherapy labeling claims are substantially equivalent to Ultrafree Max Latex Powder-Free Surgical Gloves in that they provide the following characteristics:

- intended use
- size, design, packaging
- made of natural rubber latex
- physical properties
- powder-free

Summary of Testing:

<u>Test</u>	Result
Systemic Toxicity	Glove does not elicit any toxic reactions to acute application.
Intracutaneous Reactivity	No reactivity was observed.
Guinea Pig Maximization	Glove does not display any potential for irritation.
Ultimate Elongation & Tensile Strength	Glove meets or exceeds requirements for rubber surgical gloves per ASTM D3577-99.
Barrier Defects	Glove exceeds requirements per 21 CFR §800.20 and ASTM D3577-99.
Data/Test Method	Glove meets powder level requirements for "Powder Free" designation per ASTM D3577-99.



JUN - 7 2001

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Erica Sethi Manager of Regulatory Affairs Allegiance Healthcare Corporation 1500 Waukegan Road, Building WM McGraw Park, Illinois 60085

Re: K010862

Trade/Device Name: Ultrafree Max Sterile Latex Powder-Free Surgical Gloves with Protein Labeling Claim, 50 Micrograms or Less of Total Water Extractable Protein

Per Gram Plus Chemotherapy Claim

Regulation Number: 878.4460

Regulatory Class: I Product Code: KGO Dated: March 5, 2001 Received: March 22, 2001

Dear Ms. Sethi:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the <u>Code of Federal Regulations</u>, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to

comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

Timothy A. Ulatowski

Director

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and

Radiological Health



(Division Sign-Off)

1510(k) Number ...

Division of Dental, Infection Control, and General Hospital Devices

Allegiance Healthcare Corporation 1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2460

Page 1 of 1

Applicant:	Allegiance Healthcare Corporation
510(k) Number (if kr	nown): K010862
Device Name:	Ultrafree Max Sterile Latex Powder-Free Surgical Gloves & Prefein labeling
protect a surgical wo	These gloves are intended to be worn by operating room personnel to und from contamination. These gloves contain 50 micrograms or less of total stein per gram. These gloves have been tested for use with chemotherapy
·	
(PLEASE DO NOT V	WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
Co	oncurrence of CDRH, Office of Device Evaluation (ODE)
Prescription Use(Per 21 CFR 801.109	Over-The Counter Use
Din 5. 1.	

6