

K010908 (10F4)

AUG 1 2001



**medicon**  
**Instrumente**

**Medicon eG**

**Traditional 510(k) - HCH**

MEDICON eG · POSTFACH 4455 · D-78509 TUTTLINGEN

**CHIRURGIE / DENTAL INSTRUMENTE**  
**SURGICAL / DENTAL INSTRUMENTS**

MEDICON eG, Chirurgiemechaniker-Genossenschaft  
Gänsäcker 15, D-78532 TUTTLINGEN (parcel post)  
Postfach 44 55, D-78509 TUTTLINGEN (letter post)  
GERMANY

Telefon: (0 74 62) 20 09-0  
Cables: medicon-tuttlingen

Telefax: (0 74 62) 20 09-50

e-mail: [marketing@medicon.de](mailto:marketing@medicon.de)  
[sales@medicon.de](mailto:sales@medicon.de)

homepage: <http://www.medicon.de>

Ihr Zeichen und Datum  
Your Ref. and Date

Unsere Zeichen  
Our Ref.

Telefon-Durchwahl:  
Direct dialling (phone):  
(0 74 62) 20 09-

Datum  
Date

2.

**510(k) SUMMMARY**  
**of Safety and Effectiveness**

+

**SE Comparison Table**

**Medicon eG**

[As required by Section 807.92)]

2.1

**Submitter:** [807.92 (a)(1)]

Medicon eG

Gänsäcker 15

D-78532 Tuttlingen

Germany

Tel. +49 7462 2009 0

Fax +49 7462 2009 50

eMail [sales@medicon.de](mailto:sales@medicon.de)  
[marketing@medicon.de](mailto:marketing@medicon.de)

2.2

**Date Summary Prepared:** [807.92 (a)(1)]

March 3, 2001

2-1



**2.3 Device Names:** [807.92(a)(2)]

**Proprietary** MEDICON Titanium Yasargil Aneurysm Clips

**Common** Yasargil Aneurysm Clips

**Classification** Clip, Aneurysm

**2.4. Reason for Submission:** [807.81(2)]

New Device

**2.5 Predicate Device** [807.92(a)(3)]

**Manufacturer** Aesculap Inc.  
K983758

**Proprietary Name** Yasargil Aneurysm Clip

**Catalog #'s** FT XXX T

**Manufacturer** Mizuho America Inc.  
K990202

**Proprietary Name** Sugita Aneurysm Clips

**Catalog #'s** 07 XXX-XX

**2.6 Device Description** [807.92(a)(4)]

The MEDICON Titanium Aneurysm Clips will be available as temporary or permanent devices.

The shapes, sizes, and materials used are **substantially equivalent** to those of SE devices.

**2.7 Intended Use:** [807.92 (a)(5)]

The intended use of the Medicon Titanium Yasargil Aneurysm Clips is to occlude cerebral aneurysms in either temporary or permanent manner. They are applied by MEDICON Clip Appliers.

The Aneurysm Clips are intended for use and handling by professional and qualified surgeons.

## 2.8 Environment of Use

The MEDICON Titanium Yasargil Aneurysm Clips are intended for use in healthcare facilities, including hospitals, medical clinics and surgical centers.

Within these facilities the MEDICON Titanium Yasargil Aneurysm Clips may be located in areas where sterile surgical/dental instruments are used such as operating rooms for surgery.

## 2.9 Difference in Design and Technological Characteristics when Compared to SE Devices [807.92(a)(6)]

**Material:** Like the SE devices, the MEDICON Titanium Yasargil Aneurysm Clips are manufactured from nonmagnetic, high quality Titanium Alloy (TiAl6V4) after DIN ISO 5832-3 and ISO 9713.

**Design:** The design is very similar. Both MEDICON and SE devices consist of various sizes regarding blade length and fenestration diameter and maximum opening of the jaws.

## 2.10 Discussion of Safety and Effectiveness: [807.92 (b)]

Clinical results of Aneurysm Clips made out of Titanium (TiAl6V4) and manufactured after ISO 9713 have shown their safety and effectiveness.

## 2.11 Industry Standards: [807.92 (d)]

MEDICON certifies compliance with required ISO/EN/ASTM and other device-related standards that apply to the manufacture, packaging, labeling, sterilization, and reprocessing of subject devices including the validation of these processes.

### **Substantially Equivalent**

## 2.12 Information Bearing on the Safety and Effectiveness: [807.92 (b)(3)]

The MEDICON Titanium Yasargil Aneurysm Clips have the same intended use as predicate devices. They are made of identical material. The slight differences in design and size do not adversely affect the safety and effectiveness of this Aneurysm Clips.

**The results of design validation and clinical testing raise no new issues of safety and effectiveness.**

**2.13 Comparison with predicate devices (table)**

	Medicon	Aesculap	Mizuho
<b><u>INTENDED USE</u></b>			
The intended use of Aneurysm Clips is to occlude cerebral aneurysms in either temporary or permanent manner	YES	YES	YES
<b><u>METHODS OF STERILIZATION</u></b>			
Steam Sterilization Processes (DIN 58953-9) : 134° C, 2 bar, induction time at least 5 minutes, or 121° C, 1 bar, induction time at least 15 minutes.	YES	YES	YES
<b><u>MATERIAL</u></b>			
TiAL6V4 DIN ISO 5832-3	YES	YES	YES
<b><u>DESIGN</u></b>			
Colour, Coding	YES	YES	YES
Atraumatic Jaws	YES	YES	YES
Labeling (Packaging shows indications for clip identification)	YES	YES	Assumed
Manufactured after ISO 9713 (Neurosurgical implants: Self-closing intracranial aneurysm clips)	YES	YES	Assumed



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

AUG 1 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Joachim Schmid  
General Manager  
Medicon, E.G.  
Gaensaecker 15  
D-78532 Tuttlingen,  
Germany

Re: K010908

Trade/Device Name: Medicon Titanium Yasargil Aneurysm Clips  
Regulation Number: 882.5200  
Regulatory Class: II  
Product Code: HCH  
Dated: June 6, 2001  
Received: June 13, 2001

Dear Mr. Schmid:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a stylized initial "C" and a long horizontal flourish extending to the right.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and  
Radiological Health

Enclosure

510(k) Number K010908  
Device Name **Medicon Titanium Yasargil Aneurysm Clips**  
Classification **Clip, Aneurysm**  
Product Code **84 HCH Class II 21 CFR 882.5200**

## INDICATIONS FOR USE

The MEDICON Titanium Aneurysmen Clips are to occlude cerebral aneurysms in either a temporary or permanent manner. They are applied by MEDICON Clip Appliers.

The Aneurysm Clips are intended for use and handling by professional and qualified surgeons.

---


(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)

  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K010908