

MAY - 4 2001

K011098

Special 510(k) Premarket Notification
510(k) Summary of Substantial Equivalence

Medtronic Model 9526 Reveal Plus ILR System

510(K) SUMMARY OF SUBSTANTIAL EQUIVALENCE

Proprietary Name: Reveal® Plus Insertable Loop Recorder (ILR) System.
This system is composed of the Model 9526 implantable recorder and the Model 6191 Activator. The Model 9809 Reveal Plus ILR Software, Model 9790/C programmer and programming head are also part of the system.

Common Name: Insertable Loop Recorder

Device Classification: Class II

Product Classification and Code: Cardiac Implantable Event Recorder (Product Code 74 MXC) (21 CFR 870.2800)

Contact Person: Stacey Paetschow Wessman
Product Regulation Manager
Medtronic, Inc.
4000 Lexington Avenue North
Mail Stop: X255
Shoreview, MN 55126-2983
Telephone: 763.514.9840
Facsimile: 763.514.9954
E-mail: stacey.p.wessman@medtronic.com

Predicate Device

The following device is the predicate device for the Reveal Plus ILR:

- Model 9526 Reveal Plus ILR (K994331/K003667)

Performance Standard

Performance standards do not currently exist for these devices. None established under Section 514.

Device Description

The Medtronic Reveal Plus ILR system is designed to record and store electrocardiogram (ECG) during symptomatic events. The system consists of the Model 9526 ILR and the

Model 6191 Patient Activator. A Medtronic Model 9790/C programmer equipped with a Medtronic Model 9766A/ 9766AL/9767/9767L radio frequency programming head and Model 9809 software is required for programming and retrieving data from the recorder.

Indications For Use

The Reveal Plus ILR is an implantable patient- and automatically-activated monitoring system that records subcutaneous ECG and is indicated for

- Patients with clinical syndromes or situations at increased risk of cardiac arrhythmias.
- Patients who experience transient symptoms that may suggest a cardiac arrhythmia.

Substantially Equivalent Devices

The Medtronic Reveal Plus Model 9526 ILR is believed to be substantially equivalent to the following predicate device currently in interstate commerce with respect to comparable features, materials of construction and intended use.

- Model 9526 Reveal Plus ILR (K994331/K003667)

Summary of Studies

Testing was performed in support of the this labeling modification is currently under review as part of Real-Time PMA-S for Medtronic Model 9809 Series 2.0 Reveal Plus ILR Software (Document Control Number P890003/SX, Sent to FDA on 09 April 2001). There are **no** modifications to the electrical or mechanical platform of the Reveal Plus Model 9526 ILR or the Model 6191 Patient Activator.

Conclusion (Statement of Equivalence)

Through data and information presented, numerous similarities support a determination of substantial equivalence, and therefore market clearance of the Medtronic Reveal Plus Model 9526 Insertable Loop Recorder and Model 6191 Patient Activator is supported through this Special 510(k) Premarket Notification.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY - 4 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Stacey Paetschow Wessman
Product Regulation Manager
Medtronic, Inc.
Cardiac Rhythm Management
7000 Central Avenue NE
Minneapolis, MN 55432-3576

Re: K011098

Trade Name: Model 9809 Series 2.0 Software for the Model 9526 Reveal Plus
Insertable Loop Recorder
Regulatory Class: II (two)
Product Code: 74 MXC
Dated: April 10, 2001
Received: April 14, 2001

Dear Ms. Wessman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the


Page 2 - Ms. Stacey Paetschow Wessman

Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4645. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "J. E. Dillard III", is written over the typed name.

James E. Dillard III

Director

Division of Cardiovascular

and Respiratory Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosures

INDICATIONS FOR USE

510(k) Number (if known): N/A K011098

Device Name: Reveal® Plus Insertable Loop Recorder

Indications For Use:

The Reveal Plus ILR is an implantable patient- and automatically-activated monitoring system that records subcutaneous ECG and is indicated for

- Patients with clinical syndromes or situations at increased risk of cardiac arrhythmias.
- Patients who experience transient symptoms that may suggest a cardiac arrhythmia.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

[Signature]
Division of Cardiovascular & Respiratory Devices
510(k) Number K011098

(Optional Format 1-2-96)