

Allegiance Healthcare Corporation 1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2461

SMDA REQUIREMENTS

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS Airlife® HEPA Filter

Manufacturer:

Allegiance Healthcare Corporation

1660 Iowa Avenue Riverside, CA 92507

Regulatory Affairs Contact:

Sharon Robbins

1500 Waukegan Road MPWM

McGaw Park, IL 60085

Telephone:

(847) 785-3311

Date Summary Prepared:

April, 2001

Common Name:

Airlife ® HEPA Filter

Classification:

Class II per 21CFR § 868. 5260

Predicate Device:

American Two Way Bacterial/Viral Retentive

Filter

Description:

The Airlife Breathing Circuit Bacterial Filter is a bidirectional breathing device used to reduce the

transmission of microorganisms in gases delivered to and exhaled from patients and breathing systems. The device contains a

pleated filter media that is polypropylene based media. The filter media is housed in a clear

shell.



Allegiance Healthcare Corporation 1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2461

SMDA REQUIREMENTS (continued)

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS Airlife®HEPA Filter

Intended Use:

Breathing Circuit Bacterial Filter is a device intended

to remove microbiological and PM from gases in the

breathing circuit.

Substantial Equivalence:

The Airlife®Breathing Circuit Bacterial Filter is substantially equivalent to the American Two Way Bacterial/Viral Retentive Filter in that:

- the intended use is the same

- the performance attributes are

similar

Summary of testing:

All materials used in the fabrication of the Airlife® HEPA Filter were evaluated through biological qualification safety tests as outlined in ISO 10993 Part-1 "Biological Evaluation of Medical Devices". These materials also were tested in accordance with industry recognized test methods and were found to

be acceptable for the intended use.



JUN - 8 2001

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Sharon Robbins
Allegiance Healthcare Corporation
1500 Waukegan Road
McGaw Park, IL 60085-6787

Re:

K011132

Airlife® HEPA Filter

Regulation Number: 868.5260 Regulatory Class: II (two) Product Code: 73 CAH Dated: May 18, 2001 Received: May 25, 2001

Dear Ms. Robbins:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish

Page 2 – Ms. Sharon Robbins

further announcements concerning your device in the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

James E. Dillard III

Director

Division of Cardiovascular and

Respiratory Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure



Allegiance Healthcare Corporation 1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2461

Page 1 of 1

510(k) Number (if known):	K011132	<u> </u>
Device Name:	Airlife ® HEPA FIlter	
Indications For Use:	<u> </u>	rial Filter is a device intended cal and particulate matter from circuit.
PLEASE DO NOT WRITE BEI	LOW THIS LINE - CONTIN	IUE ON ANOTHER PAGE)
Concurrence	of CDRH, Office of Device Division of Cardiovascular & Res 510(k) Number KOLL	
Prescription Use Per 21 CFR 801.109)	,,	Counter Use