

XIII. SUMMARY OF SAFETY AND EFFECTIVENESS



510(k) SUMMARY OF SAFETY AND EFFECTIVENESS POWDER-FREE STERILE SYNTHETIC EXAMINATION GLOVES

Applicant/Sponsor: Allegiance Healthcare Corporation

1500 Waukegan Road McGaw Park, IL 60085

Regulatory Affairs Contact: Erica Sethi

Allegiance Healthcare Corporation 1500 Waukegan Road, Bldg. WM

McGaw Park, IL 60085

Telephone: (847) 785-3337

Date Summary Prepared: April 12, 2001

Product Description: Powder-Free Sterile Synthetic Examination Gloves

Common Name: Examination Glove

Classification: Patient Examination Glove

Predicate Devices: Triflex Sterile Synthetic Examination Gloves

Description: These Examination Gloves are formulated using vinyl and offered sterile and powder-free.

Intended Use: These examination gloves are disposable devices intended for medical purposes that are worn on the examiner's hands or fingers to prevent contamination between patient and examiner.

Substantial Equivalence: Powder-Free Sterile Synthetic Examination Gloves are substantially equivalent to Triflex Sterile Synthetic Examination Gloves in that they provide the following characteristics:

- same intended use
- same sizes
- both made of vinyl
- both offered sterile and beaded

Summary of Testing:

| <u>Test</u> | | Result |
|-------------|--|--------|
| | | |

Primary Skin Irritation Gloves show no reactivity.

Guinea Pig Maximization Gloves do not display any potential for irritation.

Tensile Strength Gloves meet or exceed requirements per ASTM D5250-00.

Barrier Defects Gloves meet or exceed requirements per 21 CFR§800.20

and ASTM D5250-00.



MAY 2 5 2001

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Erica Sethi
Manger of Regulatory Affairs
Allegiance Healthcare Corporation
1500 Waukegan Road
McGaw Park, Illinois 60085

Re: K011186

Trade/Device Name: Powder-Free Sterile Sythetic Vinyl

Examination Gloves

Regulation Number: 880.6250

Regulatory Class: I Product Code: LYZ Dated: April 12, 2001 Received: April 18, 2001

Dear Ms. Sethi:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

/ Timoth

Timothy A. Ulatowski Director

Sugar Runner

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health Allegiance Healthcare Corporation 1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2460

Page 1 of 1

| Applicant: | Allegiance Healthcare Corporation |
|----------------|---|
| 510(k) Number: | K011186 |
| Device Name: | Powder-Free Sterile Synthetic Examination Gloves (Viny) |
| | A patient examination glove is a disposable device intended for medical on the examiner's hand or finger to prevent contamination between patient |

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

or

Over-The Counter Use X

(Division Sign-Off)

Division of Dental, Infection Control,

and General Hospital Devices 510(k) Number