

SEP 21 2001

K011341

Page 1 of 1

**510(k) SUMMARY** as required by 807.92 (c)

**APPLICANT**

**NAME:** LH Medical Products, Inc.  
**ADDRESS:** 301 E. Arrow Hwy, Suite 104  
San Dimas, CA 91773  
**PHONE:** 877-592-0523  
**FAX:** 909-592-8596  
**E-MAIL:** [lhoffstetter@lhmedicalproducts.com](mailto:lhoffstetter@lhmedicalproducts.com)

**REGISTRATION #:** 9041241

**ACTIVITY OF APPLICANT:** Initial Distributor

**CONTACT PERSON:** Lonnie Hutson  
Vice President  
**PHONE:** 877-592-0523  
**FAX:** 909-592-8596  
**E-MAIL:** [lhutson@lhmedicalproducts.com](mailto:lhutson@lhmedicalproducts.com)

**NAME OF DEVICE:**  
Trade Name: LH400 YES-SET Y-Port Enteral Spike  
Pump Set  
Common Name: Enteral Pump Set  
Classification Name: Tubes, Gastrointestinal (And  
Accessories)

**PRODUCT CODE:** KNT

**MANUFACTURER:** Weihai Medical Polymer (Group) Co., LTD  
No. 35 Yantai (W) Rd., Weihai,  
Shandong 264209  
Peoples Republic of China

**PREDICATE DEVICES:** Novartis 199307 Enteral Delivery Pump  
Set and Corpak Enteral Extension Y-  
Adapter



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SEP 21 2001

Mr. Lonnie Hutson  
Vice President  
LH Medical Products, Inc.  
301 E. Arrow Hwy., Suite 104  
SAN DIMAS CA 91773

Re: K011341  
Trade/Device Name: LH400 YES-SET Y-Port  
Enteral Spike Pump Set  
Regulation Number: 21 CFR §876.5980  
Regulation Name: Gastrointestinal tube and  
accessories  
Regulatory Class: II  
Product Code: 78 KNT  
Dated: August 18, 2001  
Received: August 21, 2001

Dear Mr. Hutson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

510(k) Number (if known): K011341

Device Name: LH400 YES-SET Y-Port Enteral Spike Pump Set

**Indications For Use:**

This device is intended to deliver, through an enteral feeding pump, liquid nutrition formulas from a container to an enteral access device (a feeding tube). It also contains a "Y" access device for enteral irrigation and medication administration without disconnecting the feeding set from the feeding tube.

(PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓ OR Over-The-Counter Use \_\_\_\_\_

Nancy C. Brogdon  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K011341