

JUN 13 2001

510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K011463

- 1. Submitter name, address, contact** Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, New York 14626-5101
(716) 453-4253

Contact Person: Darlene J. Phillips

- 2. Preparation date** Date Special 510(k) prepared: May 11, 2001
-

- 3. Device name** Trade or Proprietary Name: VITROS Immunodiagnostic Products Folate Range Verifiers
Common Name: Quality Control Material
Classification Name: Quality Control Material (Assayed)
-

- 4. Predicate device** The VITROS Immunodiagnostic Products Folate Range Verifiers (modified device) are substantially equivalent to the VITROS Immunodiagnostic Products Folate Range Verifiers (original device), (K990026, January 29, 1999).
-

Continued on next page

510(k) Summary, Continued

**5. Device
description**

The VITROS Immunodiagnostic System uses luminescence as the signal in the quantitative and semi-quantitative determination of selected analytes in human body fluids, commonly serum and plasma. Coated microwells are used as the solid phase separation system.

The system is comprised of three main elements:

1. The VITROS Immunodiagnostic Products range of immunoassay products (which in this case include VITROS Immunodiagnostic Products Folate Reagent Packs 1 and 2, VITROS Immunodiagnostic Products Vitamin B12/Folate Reagent Pack 3, VITROS Immunodiagnostic Products Red Cell Folate Reagent and VITROS Immunodiagnostic Products Folate Calibrators (cleared for market by a separate 510(k) pre-market notification K001266), VITROS Immunodiagnostic Products Anemia Controls (cleared for market by a separate 510(k) pre-market notification K 990016) and VITROS Immunodiagnostic Products Folate Range Verifiers which are combined with the VITROS Immunodiagnostic System to perform the VITROS Folate assay).
2. The VITROS Immunodiagnostic System – instrumentation, which provides automated use of the immunoassay kits. The VITROS Immunodiagnostic System was cleared for market by a separate 510(k) pre-market notification (K962919).
3. Common reagents used by the VITROS System in each assay. The VITROS Immunodiagnostic Products Signal Reagent and VITROS Immunodiagnostic Products Universal Wash Reagent were cleared as part of the VITROS Immunodiagnostic Products Total T3 Reagent Pack and VITROS Immunodiagnostic Products Total T3 Calibrators 510(k) premarket notification (K964310).

The VITROS System and common reagents are dedicated specifically for use only with the VITROS Immunodiagnostic Products range of immunoassay products.

**6. Device
intended
use**

Assayed for use in verifying the calibration range of the VITROS Immunodiagnostic System when used for the measurement of folate in human serum, plasma (heparin) and whole blood. For *in vitro* diagnostic use.

Continued on next page

510(k) Summary, Continued

7. **Comparison to predicate device** The VITROS Immunodiagnostic Products Folate Range Verifiers (modified device) are substantially equivalent to the VITROS Immunodiagnostic Products Folate Range Verifiers (original device) which were cleared by the FDA (K990026) for IVD use.

There is no change in the fundamental scientific technology of the modified device.

Table 1

Table 1 lists the characteristics of the VITROS Folate Range Verifiers (modified device) and the VITROS Folate Range Verifiers (original device).

Device Characteristic	VITROS Folate Range Verifiers (Modified device)	VITROS Folate Range Verifiers (Original device)
Storage temperature	Store at $\leq -18^{\circ}\text{C}$.	Store at $2-8^{\circ}\text{C}$.
Intended use	No change.	For use in verifying the calibration range of the VITROS Immunodiagnostic System when used for the quantitative measurement of folate.
Matrix of Range Verifiers	No change.	A buffered base matrix spiked with analyte.
Range Verifier levels	No change.	Low and high.

8. **Conclusion** The data presented in the pre-market notification demonstrate that the VITROS Folate Range Verifiers are substantially equivalent to the cleared predicate device.

The data presented in the premarket notification provide a reasonable assurance that the VITROS Folate Range Verifiers are safe and effective for the stated intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Ms. Darlene J. Phillips
Regulatory Associates
Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, NY 14626-5101

JUN 13 2001

Re: 510(k) NUMBER: K011463
Trade/Device Name: VITROS Immunodiagnostic Products Folate Range Verifiers
Regulation Number: 862.1660
Regulatory Class: I, reserved
Product Code: JJX
Dated: May 11, 2001
Received: May 14, 2001

Dear Ms. Phillips:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.

A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Statement of Intended Use

Page 1 of 1

510(k) Number (if known): K011463

Device Name: VITROS Immunodiagnostic Products Folate Range Verifiers

Indications for Use: Assayed for use in verifying the calibration range of the VITROS Immunodiagnostic System when used for the measurement of folate in human serum, plasma (heparin) and whole blood. For *in vitro* diagnostic use.

Fred Lacy
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K011463

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)