

JUL 10 2001

510(K) Sub. - Cymun-NBR Powder Free Blue Nitrile Patient Exam Glove
Submission Date: May 14, 2001
510(K) Number: _____

[TAB #9]

ATTACHMENT #7

K011497

CYMUN NBR Powder Free Blue Nitrile Patient Examination Glove

SUMMARY OF 510(k) Submission

A. INFORMATION

1. SUBMITTER'S NAME:

TILLOTSON HEALTHCARE CORPORATION

ADDRESS:

360 Route 101
Bedford, NH 03110 U.S.A.

TELEPHONE NUMBER:

(603) 472-6600

CONTACT PERSON:

Thomas N Tillotson

DATE SUMMARY PREPARED:

April 18, 2001

2. NAME OF DEVICE

TRADE OR PROPRIETARY NAME:

Cymun-NBR Powder Free Blue Nitrile
Examination Gloves

COMMON OR USUAL NAME:

Examination Glove

CLASSIFICATION NAME:

Examination Glove

3. PREDICATE DEVICE IDENTIFICATION NAME, NUMBER

1. Dual Advantage Powderfree Nitrile
Examination Glove K905765A

4. DESCRIPTION OF DEVICE

a. HOW THE DEVICE FUNCTIONS:

Nitrile Rubber Latex films form a barrier to body fluids and bloodborne
pathogens.

b. SCIENTIFIC CONCEPTS THAT FORM THE BASIS FOR THE DEVICE:

The nitrile rubber is water tight under normal conditions of use. It's tensile
properties cause it to conform to the hand, allowing movements necessary for a
medical procedure.

c. PHYSICAL AND PERFORMANCE CHARACTERISTICS SUCH AS DESIGN, MATERIALS AND PHYSICAL PROPERTIES:

Nitrile Rubber Latex is known to create a barrier to bloodborne pathogens and
and body fluids. ASTM conforming tensile properties create a glove that is strong
and flexible. The leaching process removes traces of chemical accelerants that
may be chemically irritating. The glove is manufactured in accordance with the
requirements of ASTM D 6319-00 and ASTM D5151-92 requirements.

5. STATEMENT OF INTENDED USE, INCLUDING DESCRIPTION OF THE DISEASES OR CONDITIONS THAT THE DEVICE WILL ADDRESS

This is a disposable device, intended for medical purposes, that is worn on the examiner's hand to prevent contamination between patient and examiner.
Examination gloves made with Nitrile are suitable in situations where healthcare worker or patient allergic sensitivity to Natural Rubber Latex proteins may be a factor.
Powder free gloves eliminate issues of powder contamination.

6. EXPLANATION OF SIMILARITIES OR DIFFERENCES TO PREDICATE DEVICE

- The product has similar properties compared to the predicate product.
It is suitable for situations where a powder free, non NRL glove is desirable.
- It is blue compared to the white color of the predicate product

B. IF SE DECISION BASED ON PERFORMANCE DATA

1. DISCUSSION OF NON-CLINICAL TESTS

SPECIFICATION	PROPOSED Powder Free	PREDICATE Powder Free
	Blue Color	White Color
PERFORMANCE STANDARDS	ASTM D	ASTM D
WATER TIGHTNESS	ASTM D5151-99	ASTM D5151-99

2. DISCUSSION OF CLINICAL TESTS

SPECIFICATION	PROPOSED	PREDICATE
<u>SAFETY</u>		
RABBIT IRRITATION	Passes	Passes
GUINEA PIG SENSITIZATION	Passes	Passes

DISCUSSION OF SAFETY OR EFFECTIVENESS DATA OBTAINED

with specific reference to adverse effects and complications

3. CONCLUSIONS DRAWN FROM NONCLINICAL AND CLINICAL TESTS THAT
DEMONSTRATE
SAFETY EFFECTIVENESS, AND PERFORMANCE => PREDICATE PRODUCT

The Cymun-NBR, Examination Glove has been carefully compared to legally
marketed devices in the 510(k). The data summaries indicate that the proposed
product meets or exceeds acceptable scores for the predicate product in nonclinical
tests, and satisfies the requirements for a safe and effective powder free, nitrile
medical glove.

Pursuant to 21 C.F.R. 807.87 (j), I, Thomas N Tillotson, CEO,
certify that to the best of my knowledge and belief and based upon the data
and information submitted to me in the course of my responsibilities as the
CEO for TILLOTSON HEALTHCARE CORPORATION,
and in reliance thereupon, the data and information submitted in this
premarket notification are truthful and accurate and that no facts material to a review
of the substantial equivalence of this device have been knowingly omitted from this
submission.



Thomas N Tillotson
CEO



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 1 0 2001

Mr. Thomas N. Tillotson
Tillotson Healthcare Corporation
360 Route 101
Bedford, New Hampshire 03110

Re: K011497
Trade/Device Name: Cymun-NBR Powder-Free Blue Nitrile
Examination Gloves
Regulation Number: 21 CFR 880.6250
Regulatory Class: I
Product Code: LZA
Dated: May 14, 2001
Received: May 15, 2001

Dear Mr. Tillotson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 -Mr. Tillotson

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

- 3.0 **Indications for Use Statement:** Include the following or equivalent Indications for Use page.
The information, data and labeling claims in the entire the 510(k) submission must support and agree with the Indications for Use statement.

INDICATIONS FOR USE

Applicant: Tillotson Healthcare Corporation


5 10(k) Number (if known):* K011497

Device Name: Cymun-NBR Powder Free Blue Nitrile Examination Glove

Indications For Use:

The *Cymun-NBR*, Examination Glove is "a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner."
(21CFR 880.6250).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K011497

Prescription Use _____ OR Over-The-Counter _____
Per 21 CFR 801.109
(Optional Format 1-2-96)

For a new submission, do NOT fill in the 510(k) number blank.