

**510(k) Summary of Safety & Effectiveness**

---

**Submitter** Vanguard Medical Concepts, Inc.  
5307 Great Oak Drive  
Lakeland, FL 33815

---

**Contact** Mike Sammon, Ph.D.  
Director, Research & Development  
(863) 683-8680, extension 228  
mikes@safe-reuse.com

---

**Date** May 31, 2001

---

**Devices**

- Trade Names: Vanguard Reprocessed Electric Biopsy Forceps  
⇒ Boston Scientific Microvasive® Radial Jaw® 3, and  
⇒ C.R. Bard Precisor® Electric Biopsy Forceps
- Common Name: Electric biopsy forceps
- Classification: 21 CFR 876.4300 – Endoscopic Electrosurgical Unit and Accessories – Class II
- Product Code KGE

---

**Predicate Devices**

- Boston Scientific Microvasive® Radial Jaw® 3 Electric Biopsy Forceps believed to be legally marketed under 510(k) premarket notification K860366
- C.R. Bard Precisor® Hot Biopsy Forceps believed to be legally marketed under 510(k) premarket notification K912601

---

**Indications for Use** When used with a compatible electrosurgical unit, endoscope and patient grounding pad, electric biopsy forceps are intended for electrocautery and removal of polyps and/or tissue within the gastrointestinal tract. The forceps are prescription devices intended for a single patient use only.

---

*Continued on next page*

## 510(k) Summary of Safety & Effectiveness, Continued

---

**Device Description**

Biopsy forceps are devices used to collect tissue samples and/or remove polyps within the gastrointestinal tract via the operating channel of endoscopic instruments.

Electric biopsy forceps include an electrocautery function when used with a compatible electrosurgical unit.

The devices consist of flexible, electrically isolated sheaths with distal grasping cups that are controlled by a proximal handle. The forceps are guided by endoscopy through a biopsy channel with a minimum dimension of 2.8 mm. Monopolar electrocautery requires electrical connection of the forceps to a compatible electrosurgical unit and use of an appropriate patient grounding pad.

Vanguard receives previously used electric biopsy forceps from healthcare facilities; cleans, inspects, tests, repackages and sterilizes the forceps; and returns them to the healthcare facility.

---

**Technological Characteristics**

The technological characteristics of the Vanguard reprocessed electric biopsy forceps are the same as the predicate devices, the original equipment manufacturer (OEM) electric biopsy forceps. The materials and dimensions are unchanged and the physical characteristics, performance specifications, and all other characteristics are essentially identical.

---

**Test Data**

Decontamination and cleaning, packaging and sterilization validations and functional/performance, shelf life and biocompatibility testing demonstrates that the reprocessed devices perform as intended and are safe and effective.

---

**Conclusion**

Based on the information provided herein and the 510(k) "Substantial Equivalence" Decision Making Process Chart, we conclude that the Vanguard reprocessed electric biopsy forceps are substantially equivalent to the predicate device, the respective OEM biopsy forceps under the Federal Food, Drug and Cosmetic Act.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 27 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mike Sammon, Ph.D.  
Director, Research and Development  
Vanguard Medical Concepts, Inc.  
5307 Great Oak Drive  
LAKELAND FL 33815

Re: K011800  
Trade/Device Name: Vanguard Reprocessed  
Electric Biopsy Forceps  
Regulation Number: 21 CFR §876.4300  
Regulation Name: Endoscopic electrosurgical unit  
and accessories  
Regulatory Class: II  
Product Code: 78 KGE  
Dated: September 7, 2001  
Received: September 10, 2001

Dear Dr. Sammon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

K011800

### Indications for Use

510(k) Number: K011800

Device Name: Vanguard Reprocessed Electric Biopsy Forceps

#### Indications for Use:

When used with a compatible electrosurgical unit, endoscope and patient grounding pad, electric biopsy forceps are intended for electrocautery and removal of polyps and/or tissue within the gastrointestinal tract. The forceps are prescription devices intended for a single patient use only.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED.)

-----  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)

Nancy C Brogdon  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K011800