Medical Instruments Technology, Inc's. Reprocessed Trocar Premarket Notification



KO12659

NOV 0 8 2001

Quality Reprocessing and Surgical Cost Containment Systems

Section 12: 510k Summary

Name of Submitter

Medical Instruments Technology, Inc. 385 North 3050 East Saint George, UT 84790 Tel: (435) 674-4010

Fax: (435) 674-9819

Contact persons

Tom Haueter, RA/QA Manager Crystal Batcabe, Assistant RA/QA Manager

Summary Prepared August 10, 2001

Device Name and Classification

Common Name:

Trocar

Classification:

Class II per 21 CFR 876.1500

Predicate Device

Ethicon Trocar (K971475)

Description of Device

MIT's reprocessed Trocars are composed of two main parts, an inner piece and an outer piece. The inner piece has a sharp blade and a spring-loaded shield or a blunt edge. The blade (or blunt edge) allows for insertion into the operative cavity. The shield provides for protection of internal organs from puncture once the operative cavity has been entered. The shield is automatically engaged when the operative cavity has been entered. The outer piece consists of a cannula with a gasket. The cannula allows for introduction of surgical instruments while the gasket assures pneumoperitoneum during the procedure.

Intended Use

The device is intended to establish a path for the introduction of minimally invasive instruments. Trocars, which are provided with a stopcock, are also intended to insufflate the operating area.

Medical Instruments Technology, Inc's. Reprocessed Trocar Premarket Notification

Technological Characteristics

MIT's reprocessed trocars have the same technological characteristics as the predicate devices. MIT does not change any of the design characteristics or materials during reprocessing. MIT has shown that the reprocessed trocars are substantially equivalent to the predicate devices by performance of the physical and functional tests. Additionally, we have tested the device for biocompatibility by performing the bioburden test and the ETO residual test. In all tests, the reprocessed devices have been equal the predicate devices



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

NOV 0 8 2001

Mr. Jack Speer President Medical Instruments Technology, Inc. 385 North 3050 East Suite B St. George, Utah 84790

Re: K012659

Trade Name: Dilating Tip Trocar Regulation Number: 876.1500

Regulation Name: Laparoscope, General/Plastic Surgery

Regulatory Class: Class II

Product Code: GCJ Dated: August 10, 2001 Received: August 13, 2001

Dear Mr. Speer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Jack Speer

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Sisan walky M

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

510(k) Number (if known): <u>K 0126.57</u>	<u></u>	
Device Name: Cannulated Trocass		_
Indications For Use:	•	٠.

MIT's reprocessed trocars are indicated for use to establish a path of entry for minimally invasive instruments.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division o General, Restorative and Neurological Devices

510(k) Number <u>K01265</u> Over-The-Counter Use

OR

Prescription Use V (Per 21 CFR 801.109)