

DEC 20 2001

K012660

Pg 1 of 4

## 2. 510(k) SUMMARY of Safety and Effectiveness

### GIMMI GmbH

As required by Section 807.92(c)

#### 2.1 Submitter: [807.92 (a)(1)]

GIMMI GmbH  
Junkerstr. 1  
D-78532 Tuttlingen  
Germany

Tel. +49-7461 - 9 65 90-0  
Fax +49-7461 - 9 65 90-33  
e-mail [info@gimmi.de](mailto:info@gimmi.de)

#### 2.2 Contact Person: [807.92 (a)(1)]

Dagmar S. Mäser  
Business Support International  
Amstel 320-I  
1017 AP Amsterdam  
The Netherlands

Tel. +31-20-428 95 91  
Fax +31-20-428 94 29  
eMail [bsi@xs4all.nl](mailto:bsi@xs4all.nl)

#### 2.3 Date Summary Prepared: [807.92 (a)(1)] December 10, 2001

#### 2.4 Device Names: [807.92 (a)(2)]

**Proprietary** GIMMI ALPHA® gastro-urology and laparoscopic endoscopes, endoscopic accessories and GIMMI devices for minimally invasive GI, GU, and laparoscopic diagnostic and/or therapeutic indications.

**Common** Endoscopic Instruments & Accessories and Minimally Invasive GI and GU Devices

#### PANEL 78

Cystoscope, Diagnostic	78 FAJ	876.1500	II
Needle, Endoscopic	78 FBK	876.1500	II
Resectoscope & Accessories	78 FDC	876.1500	II
Set, Laparoscope	78 FDE	876.1500	II
Sheath, for Endoscope	78 FED	876.1500	II
Urethroscope	78 FGC	876.1500	II
Resectoscope	78 FJL	876.1500	II
Laparoscope, General & Plastic Surgery	78 GCJ	876.1500	II
Endoscopes & Accessories	78 K0G	876.1500	II
Urethrotome	78 EZO	876.4770	II
Electrode, Electrosurgical, Active, Urol	78 FAS	876.4300	II
Unit, Electrosurgical,	78 KNS	876.4300	II
<b>Minimally Invasive GI and GU Devices</b>			
Evacuator, Gastro-Urology	78 KQT	876.4370	II
Cannula and Trocar, Suprapubic	78 FBM	876.5090	II

Dislodger, Stone Biliary	78LQR	876.5010	II
Clamp, Penile	78 FHA	876.5160	II
Catheter, Urological	78 KOD	876.5130	II
Dilator, Urethral	78 KOE	876.5520	II

**PANEL 79**

Device, Electros., Cutting & Coag'n & Acc	79 GEL	878.4400	II
--	--------	----------	----

**PANEL 78: EXEMPT DEVICES**

Forceps, Biopsy, Non-Electric	78 FCL	876.1500	I Exempt
Accessories, Cleaning	78 MNL	876.1500	I Exempt
Brushes for Endoscope			
Surg Instruments, G-U, Manual & Acces	78 KOA	876.4300	I Exempt
Bougies, Urological	78 FAX	876.5520	I Exempt

**PANEL 79: EXEMPT DEVICES**

Catheter, Cholangiography	79 GBZ	878.4200	I Exempt
Tray, Surgical, Instrument	79 FSM	878.4800	I Exempt
Cannula, Surgical, General & Plastic Surg	79 GEA	878.4800	I Exempt

**2.5 Reason for Submission:**  
New Devices

**2.6 Predicate Devices: [807.92 (a)(3)]**

Predicate devices are produced by

Günter Bissinger Medizintechnik

Comeg Endoscopy

Dufner Instrumente GmbH

Henke-Sass Wolf, GmbH

Optus, Inc.

Pilling Weck Group

Wolf

Karl Storz Endoscopy

and a wide range of other manufacturers, including:

Jarit (J. Jamner Surgical Instruments, Inc.)

Snowden-Pencer, Inc.

SurgiTech, Inc. (Surgical Technologies International, Inc.)

United States Surgical Corp.

Allegiance Healthcare Corp.

**2.7 Device Description: [807.92(a)(4)+(6)]**

GIMMI ALPHA Endoscopes are comprised of rigid, panoramic telescopes using rod lens technology. The body contact portions are composed of surgical grade stainless steel, which is commonly used in medical devices for a wide range of applications and has a long history of biocompatibility for human use.

K 012 660  
Pg 3 of 4

**Laparoscopic and urological-gastroenterologic and related accessories** are composed of reusable handle and shaft assemblies and removable, reusable tip assemblies. Needle holders and other Class I devices included in the endoscopic catalogs may be one-piece. The instruments are designed and manufactured specifically for the purpose of manipulating soft tissue structures (grasping, cutting, dissecting, coagulating and suturing).

**Minimally Invasive GI and GU Devices** are design specific for short term use in diagnostic and/or therapeutic procedures

**2.8 Intended Use:** [807.92 (a)(5)]

GIMMI ALPHA® gastro-urology and laparoscopic endoscopes, endoscopic accessories and GIMMI devices for minimally invasive GI, GU, and laparoscopic diagnostic and/or therapeutic indications are intended to be used by qualified physicians to provide access, illumination and visualization of internal structures and for manipulating soft tissue (grasping, cutting, coagulating, dissecting, and suturing) in a wide variety of diagnostic and therapeutic laparoscopic/urologic closed and minimally invasive procedures.

**2.9 Industry Standards/Performance Data:** [807.92 (d)]

GIMMI certifies compliance with relevant ISO/EN/ASTM/AAMI/ANSI/IEC and other device-related standards that apply to the manufacture, packaging, labeling, and reprocessing of subject devices including the validation of these processes.

**2.10 Summary of Testing**

All materials used in the composition of GIMMI ALPHA® gastro-urology and laparoscopic endoscopes, endoscopic accessories and GIMMI devices for minimally invasive GI, GU, and laparoscopic diagnostic and/or therapeutic indications were subjected to performance and physical tests to evaluate safety, effectiveness, and reliability of the devices. All results were in conformance with the cited harmonized device standards.

**2.11 Information Bearing on the Safety and Effectiveness:**

[807.92 (b)(3)]

The GIMMI ALPHA® gastro-urology and laparoscopic endoscopes, endoscopic accessories and GIMMI devices for minimally invasive GI, GU, and laparoscopic diagnostic and/or therapeutic indications have the same intended use as predicate devices. They are made of the same material and produced to the same international and FDA-recognized standards. Slight modifications in design do not adversely affect the safety and effectiveness of these devices.

In summary, the

- intended use
- performance attributes
- materials and
- basic design

are identical (endoscopes and HF accessories) and identical/  
substantially equivalent to SE devices.

**The results of design validation raise no new issues of safety  
and effectiveness.**



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 20 2001

GIMMI GmbH  
% Ms. Dagmar S. Mäser  
FDA Liaison  
Business Support International  
1017 AP Amsterdam  
Amstel 320-I  
THE NETHERLANDS

Re: K012660

Trade/Device Name: GIMMI *ALPHA*<sup>®</sup>  
Gastro-Urology, & Laparoscopic Endoscopes,  
Endoscopic Accessories and GIMMI Devices  
for Minimally Invasive GI, GU &  
Laparoscopic Diagnostic and/or Therapeutic  
Indications

Regulation Numbers: See Enclosure II

Regulation Names: See Enclosure II

Regulatory Class: See Enclosure II

Product Code: See Enclosure II

Dated: October 17, 2001

Received: October 19, 2001

Dear Ms. Mäser:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see Attachment I) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number in Attachment I:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892. 2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
for Nancy C. Brogdon

Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosures

- (1) Indications for Use
- (2) Device listing

510(k) Number

K012660

Device Name

**ALPHA® Endoscopic Instruments  
& Accessories**

## INDICATIONS FOR USE

GIMMI ALPHA® Gastro-Urology and Laparoscopic Endoscopes, Endoscopic Accessories and GIMMI devices for minimally invasive gastrointestinal (GI), genitourinary (GU), and laparoscopic diagnostic and/or therapeutic indications are intended to be used by qualified physicians to provide access, illumination and visualization of internal structures and for manipulating soft tissue (grasping, cutting, coagulating, dissecting and suturing) in a wide variety of diagnostic and therapeutic laparoscopic/urologic and minimally invasive procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRE, Office of Device Evaluation (ODE)

Prescription Use ☒

OR

Over-The-Counter Use ☐

(Per CFR 801.109)

(Optional Format 1-2-96)

  
(Division Sign-Off)Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K012660