

AUG 1 1996

K955534

510(k) Summary  
*Segue*<sup>TM</sup> Infusion Catheter

Trade Name: *Segue*<sup>TM</sup> Infusion Catheter

Classification Name: Catheter, Percutaneous

Classification: Class II

Submitted By: Interventional Innovations Corporation  
2670 Patton Road  
St. Paul, MN 55113  
(612) 636-6634

Contact: Karen Peterson  
Director of Clinical and Regulatory Affairs

Predicate Device: FasTRACKER® Infusion Catheter  
Target Therapeutics

Device Description

The *Segue*<sup>TM</sup> Infusion Catheter is a single lumen, over-the-wire device designed to locally deliver fluids into the peripheral (non-coronary, non-cerebral) vasculature.

Intended Use

The *Segue*<sup>TM</sup> Infusion Catheter is indicated for the delivery of therapeutic agents into the peripheral vasculature. It is not indicated for use in coronary or cerebral vasculature. It is not intended for use with thrombolytics. It is not intended for use with power injection pumps.

Testing

Physical testing of the product under simulated conditions included: dimensional inspection, deployment and recoil verification, marker band attachment, infusion flow rate, internal pressurization, bond strength, flexural fatigue strength, radial force and trackability. All testing results were within product engineering and marketing specifications.

Biocompatibility testing was performed on the sterile materials used in the construction of this infusion catheter. All materials passed the biocompatibility testing and are suitable for this application.

Animal studies were conducted to assess placement of the device in a vessel as well as thrombus accumulation on the device. There were no adverse results reported.

#### Summary of Substantial Equivalence

The *Segue*™ Infusion Catheter is constructed of the same or substantially equivalent materials as found in the predicate device. The sizes and configurations available are comparable as is the packaging methods and materials. The clinical indications for use are substantially equivalent to those of the predicate device. Because of the similarities in materials, construction, indications for use, packaging and testing results, this product does not raise any new safety or effectiveness issues.