

K960885

Physiometrix, Inc.
510(k), Premarket Notification
HydroSpot System, Model 1500 & 1501

SUMMARY OF SAFETY AND EFFECTIVENESS

Date: March 1, 1996

Company: Physiometrix, Inc.
Five Billerica Park
101 Billerica Avenue
N. Billerica, MA 01862

Contact: Dawn E. Frazer
Director, Regulatory Affairs & Quality Assurance
(508) 670-2422
(800) 474-9746

Subject Device: Model 1240 , Small e-Net (42 - 48 cm)

Predicate Device: K930080, HydroDot NeuroMonitoring System

Classification: Class II, CFR 21 Part 882.1320, Cutaneous Electrodes

Description: The small e-Net is designed to fit patients with a head size in the range of 42 - 48 cm in circumference. It is constructed using the same materials and manufacturing processes as the larger sized net 54 - 62 cm) which is described in K930080.

Intended Use: The e-Net serves two functions as a part of the Hydro Dot NeuroMonitoring System. First, it locates the EEG electrodes according the 10-20 International System. Secondly, it conducts the electrical signal sensed by the electrodes from the skin to the EEG equipment.

Design/Materials: The small e-Net is fabricated from the same materials as the larger version. Only the size of the components was reduced to achieve the proportional fit on the intended head sizes.