

A.

JUN - 3 1996

K960941

**1. Submittals:**

**Base Ten Systems, Inc.**  
One Electronics Drive  
Trenton, New Jersey 08619-0151  
Voice: (609) 586-7010  
Fax: (609) 586-1593  
Contact Person: Richard J. Farrelly  
Date of Preparation: 13 May 1996

**2. Name of Device:**

Trade or Proprietary Name: uPACS™  
Company or Usual Name: Picture Archiving System  
(Established Name)  
Classification Name: Device, Digital Image Storage - 90LMB  
Digitizer, Image - 90LMA

**3. Predicate Devices:**

Intrascan  
and  
Instrascan II

Manufactured by:

Mark Solutions, Inc.  
Mark Care Medical Systems, Inc.  
87 Route 17 North  
Maywood, New Jersey 07607  
Telephone: (201) 368-8118

**4. Device Description:**

uPACS is a software controlled microprocessor based archiving system using installed proprietary software for recording, storage and presentation of video ultrasound images. The system is designed for connection to and use with FDA cleared ultrasound imaging equipment which provides video output and is cleared for operation and use in the patient environment.

In the uPACS standard configuration, an FDA approved ultrasound imaging device (scanner) is connected to an IBM-compatible personal computer (PC). Images captured by

---

the scanner are displayed on the monitor of the PC along with any text that is typed on the PC keyboard. The ultrasound image and patient text data can be stored on a CD-ROM in the PC, and images can be exported to an FDA approved printer by the PC.

uPACS enables users to perform the following system functions:

- ❖ View live images from user provided FDA approved ultrasound scanners or other NTSC or PAL video sources;
- ❖ Capture single-frame images from ultrasound scanners or other NTSC or PAL video sources;
- ❖ Obtain text data directly from captured images;
- ❖ Edit text and image data before permanent storage;
- ❖ Review stored images on the system and on CD-ROM;
- ❖ Display up to four images at a time for review;
- ❖ Place notations directly on images;
- ❖ Archive images using CD-ROM;
- ❖ Export images to diskette and user provided printer(s);
- ❖ Provide a means for user training and proficiency demonstration;
- ❖ Assign security privileges to specific users in accordance with clinical access criteria.

**5. *Intended Use Statement:***

uPACST<sup>TM</sup> is a software controlled microprocessor based archiving system using installed proprietary software for recording, storage and presentation of video ultrasound images. The system is designed for use with FDA approved ultrasound imaging equipment which provides video output and is approved for operation and use in the patient environment.

**6.** The technological characteristics of the uPACST<sup>TM</sup> Device are essentially the same as that of the predicate device.

**B.** The determination of substantial equivalence is not based on an assessment of performance data.