

JUN 14 1996

K961481

## 510(k) Summary

- Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.
- 
- 1) Submitter name, address, contact** Boehringer Mannheim Corporation  
9115 Hague Rd.  
Indianapolis, IN 46250  
(317) 845-2327
- Contact Person: John D. Stevens
- Date Prepared: April 17, 1996
- 
- 2) Device name** Proprietary name: Elecsys 2010<sup>®</sup> analyzer
- Common name: immunoassay analyzer
- Classification name: Discrete photometric analyzer for clinical use
- 
- 3) Predicate device** We claim substantial equivalence to the Boehringer Mannheim ES300 Immunoassay system.
- 
- 4) Device Description** The Elecsys 2010 analyzer is a fully automated, random access, computer controlled analytical system for quantitative and qualitative determinations of analytes in body fluids
- 
- 5) Intended use** The Elecsys 2010 analyzer is intended to be used for the in vitro quantitative and qualitative analysis of analytes in body fluids.
- 

*Continued on next page*

## 510(k) Summary Elecsys® 2010 Analyzer, Continued

### 6) Comparison to predicate device

The Boehringer Mannheim Elecsys 2010 analyzer is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed ES300 Immunoassay analyzer.

The following table compares the Elecsys 2010 with the predicate device. Specific data on the performance of the system have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device is provided in attachment 6.

#### Similarities:

Feature	Elecsys 2010 Analyzer	ES300 Analyzer
Intended use	Intended for the in vitro quantitative/qualitative determination of analytes in body fluids	Same
Number of reagents	Up to 2	Same
Sample detection	Liquid level detection and clot detection	Same
Sample type	Serum and plasma	Same
Host interface	RS232C bidirectional	Same

#### Differences:

Feature	Elecsys 2010 analyzer	ES300 analyzer
Operating principle	Electrochemiluminescence immunoassay method type	Enzyme-linked immunosorbent assay method type
Detection system	Flow through electrochemiluminescence detection cell	Flow through cuvette with halogen lamp light source
Sample positions	30 positions for samples, controls and calibrators	150 positions for samples, controls and calibrators
Sample volume / test	10 to 50µL	5 to 200 µL
Reagent positions	15 reagent positions 2 diluent/reagent positions 1 BlankCell position	12 reagent positions 1 cleaning solution 2 universal substrate

*Continued on next page*

## 510(k) Summary, Elecsys 2010® Analyzer Continued

6) Comparison to predicate device, (cont.) Performance Characteristics:

Feature	Elecsys 2010 vs. ES300 Least Squares Regression
Method Comparison T4	$Y = 1.02x + 0.85$ SEE = 5.39 $r = 0.947$
Method Comparison FT4	$Y = 0.954X + 0.18$ SEE = 1.11 $r = 0.981$
Method Comparison TSH	$Y = 1.09X + 0.14$ SEE = 0.798 $r = 0.991; n = 132$
Method Comparison T-Uptake	$Y = 0.99X - 0.03$ SEE = 0.04 $r = 0.908; n = 319$
Method Comparison hCG	$Y = 1.35X - 9.21$ SEE = 17.50 $r = 0.989; n = 64$
Method Comparison Troponin-T	$Y = -0.033 + 0.996$ SEE = 0.547 $r = 0.969; n = 54$