



AUG 23 1996

P.O. Box 708  
Warsaw, IN 46581-0708  
219 267-6131

Summary of Safety and Effectiveness

K962196

The summary of safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990.

- Submitted by:

Zimmer, Inc.  
P.O. Box 708  
Warsaw, IN 46581-0708

- Prepared by:

Ruth Ann Wood  
Senior Regulatory Affairs Associate

- Date:

June 5, 1996

- Trade Name:

*NexGen®* Complete Knee Solution with *Co-Nidium™* Surface Hardening Process  
and *Legacy™* Knee System with *Co-Nidium™* Surface Hardening Process

- Classification:

Class II  
21 CFR 888.3530

- Predicate Devices:

Cobalt-Chromium Molybdenum Alloy Femoral Heads with *Co-Nidium* Surface Hardening Process; *MG II®* Total Knee System with *Ti-Nidium®* Surface Hardening Process; *MG II®* Porous Total Knee System, *NexGen®* Complete Knee Solution; and the *Legacy™* Knee System

- Device Description:

The femorals are made of cast *Zimaloy®* Cobalt-Chromium-Molybdenum Alloy. They are hardened by the *Co-Nidium™* Surface Hardening Process described in



They are hardened by the *Co-Nidium™* Surface Hardening Process described in Section XI. They are available in left and right configurations. *NexGen/Legacy* Knee femorals cover a multitude of patient specific needs with the Cruciate Retaining (CR), Posterior Stabilizing (PS and L-PS), Cruciate Retaining Augmentable (CRA), and Constrained Condylar Knee (L-CCK) designs. Certain femorals can be augmented with stem extensions. The extensions come in a broad variety of straight and offset shapes as well as numerous lengths and diameters. These femorals will be available in precoat, porous, and option (NoCoat) versions.

- Intended Use:

Knee femoral components are single use devices implanted in the human knee during total knee arthroplasty. This knee is intended for cement use only.

- Performance Data:

The performance data submitted in the 510(k) demonstrated the improved wear resistance compared to nontreated femoral components.

RA05603K.510