

Appendix A. 510(k) Summary of Safety and Effectiveness**510(k) Summary of Safety and Effectiveness**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: _____

Applicant Information:

Date Prepared: June 26, 1996

Name: Heartport, Inc.
Address: 200 Chesapeake Drive
Redwood City, CA 94063

Contact Person: Robert J. Chin
Phone Number: (415) 306-7900
Fax Number: (415) 306-7905

Device Information:

Trade Name: Endoaortic Clamp (Catalog Number 02055)
Common Name: Endoaortic Clamp Catheter
Classification Name: Cardiopulmonary bypass vascular catheter

Equivalent Devices:

Heartport Endoaortic Clamp (Catalog Number 01055)

Intended Use:

Occlusion of the aorta, delivery of cardioplegic solution, and monitoring of aortic root pressure during cardiopulmonary bypass.

Comparison To Predicate Devices:

The modifications to the predicate device provides the surgeon with an endoaortic clamp already preassembled with the necessary accessories in one sterile package.

Non-clinical Test Results:

Performance testing has demonstrated with 95% confidence that the Endoaortic Clamp will meet or exceed Heartport, Inc. performance standards.

Test Conclusions:

Performance testing has demonstrated that the Endoaortic Clamp will function safely and efficaciously, while meeting the anticipated clinical requirements for the intended use.