

Heelbo, Inc.

1134 N. Homan Ave.
Chicago, IL 60651
312/489-7711

Telex RCA 297-248
800/323-5444
Fax 312/489-7744

K962839

OCT 18 1996

Heelbo, Inc.
Heelbo Quilted Limb Holder with Buckle and Strap

Safety and Effectiveness Summary

1. Submitter's name, Address and Contact Person

Submitter

Heelbo, Inc.
1134 N. Homan Ave.
Chicago, IL 60651

Contact Person

Joseph S. Tokarz; Manager, Regulatory Affairs
Hollister Incorporated
2000 Hollister Drive
Libertyville, IL 60048
Ph: (847)680-2849, Fax: (847)918-3860

Date Summary Prepared - July 17, 1996

2. Name of Device:

Heelbo Quilted Limb Holder with Buckle and Strap

3. Name of Predicate Device(s)

Heelbo Quilted Limb Holder with Buckle and Strap

4. Description of Device

The Heelbo Quilted Limb Holder with Buckle and Strap is a quilted limb holder with a Velcro® fastener, buckle and straps that secure to a bed, gurney, chair or wheelchair. The Heelbo Quilted Limb Holder with Buckle and Strap is intended to be a protective limb restraint that limits the patient's movement to the extent necessary for treatment, examination, or protection of the patient or others. The adjustable buckle allows for supine and sitting positions. The Heelbo Quilted Limb Holder is machine washable.

The Heelbo Quilted Limb Holder with Buckle and Strap is available in the following sizes:

Stock No. 9300

Small-Regular / cuff size: 11 1/2" long
strap: 3/4" wide, 60" long

Stock No. 9310

Large-XL / cuff size: 13 1/2" long
strap: 3/4" wide, 60" long

5. Statement of Intended Use

Heelbo Quilted Limb Holder with Buckle and Strap is intended to be a protective limb restraint that limits the patient's movement to the extent necessary for treatment, examination, or protection of the patient or others. The device can be used in the supine or sitting positions.

Heelbo, Inc.
Heelbo Quilted Limb Holder with Buckle and Strap

6. Statement of Technological Characteristics of the Device

Heelbo Quilted Limb Holder with Buckle and Strap is a quilted limb holder with a Velcro® fastener, buckle and straps that secure to a bed, gurney, chair or wheelchair. The Heelbo Quilted Limb Holder with Buckle and Strap is intended to be a protective limb restraint that limits the patient's movement to the extent necessary for treatment, examination, or protection of the patient or others. The adjustable buckle allows for supine and sitting positions. The Heelbo Quilted Limb Holder is machine washable.

The subject devices are identical in intended use, design, materials, manufacturing process, physical and mechanical specifications and issues of safety and effectiveness to the devices prior to the submission of this notification. The only difference is that the product labeling has been revised to comply with the Agency's labeling requirements set forth in the draft "Guidance on the Content of Premarket Notification [501(k)] Submissions for Protective Restraints" dated December 1995.

7. Biocompatibility Assessment

The subject devices are identical in component materials to the predicate devices. The suppliers of the materials use in the fabrication of these devices have stated that there is a history of safe use of their materials in the clothing and garment industry. In addition, Heelbo, Inc., is not aware of any reports or complaints of skin irritation associated with these materials.

8. Conclusion

Based upon the information presented above it is concluded that the proposed Heelbo Quilted Limb Holder with Buckle and Strap are safe and effective for their intended use and are substantially equivalent to the predicate device.