Summary of Safety and Effectiveness

Owner:

Gallini s.r.l.

88. Via S. Faustino 41037-Mirandola -Italy

Contact Person:

Caterina E. Galluzzo President, Gallini Intl.

Submission Date: July 1996

Device Name:

TX® Biopsy Needle

Predicate Device: Promedical Ltd.. - PRO-B Biopsy Needle (510 (k) # K951598)

Description of

This device is an automatic disposable biopsy Needle developed for taking a range of tissue specimens, including those from the prostate,

Device function:

kidney and liver.

Intended Use:

This device is an automated, disposable and adjustable biopsy instrument designed for ultimate accuracy in obtaining high quality histological core samples, with three different choice of penetration depth and obtaining

histological samples from 0.8 to 1.9 mm length.

COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

The TX® Biopsy Needle is equivalent to the predicate devices in term of technological characteristics

The following tests was performed to show the safety, efficacy and performance of the product.

TEST Gun Firing Test (TX/F/T) see attch. protocol Pyrogen Test (LAL)	RESULT passed passed
All biocompatibility tests on the components are performed by the supplier of the row material.	